Acumed®

ACUTRAK® , ACUTRAK2® SCREW SYSTEMS

OBSOLETE
Visit www.acumed.net/ru
for the latest version.

MediMark® Europe Sarl.

11 rue Emile ZOLA. BP 2332
38033 GRENOBLE CEDEX 2
FRANCE
+33 4 76 86 43 22

EC REP

Acumed® LLC

5885 NW Cornelius Pass Road
Hillsboro, OR 97124-9432
+1 (503) 627-9957
www.acumed.net

PKGI-27-F
EFFECTIVE 4/2015

STERILE R

REF

LOT

0473

* P K G I - 2 7 - F *

TABLE OF CONTENTS

Click on language to navigate to page.

English – US	3
Dansk – DA	9
English – EN	16
Deutsch – DE	23
Ελληνικά – EL	31
Español – ES	39
Français – FR	47
Italiano – IT	55
Nederlands – NL	62
Norsk – NO	69
Português – PT	76
Suomi – FI	85
Svenska – SV	90
Türkçe – TR	97

OBSOLETE
Visit www.acumed.net/ifu
for the latest version.

US ACUTRAK®, ACUTRAK 2® SCREW SYSTEMS

FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

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Although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be conveyed to the patient.

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WARNINGS: For safe effective use of the implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for the device. The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Device breakage or damage can occur when the implant is subjected to increased loading associated with delayed union, nonunion, or incomplete healing. Improper insertion of the device during implantation can increase

the possibility of loosening, or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant. These cautions include the possibility of the device or treatment failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing, and the possibility of nerve or soft tissue damage related to either surgical trauma or the presence of the implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail. The components of these systems have not been tested for safety, heating, or migration in the MRI environment. The implants may cause distortion and/or block the view of anatomic structures on radiographic images. Similar products have been tested and described in terms of how they may be safely used in post-operative clinical evaluation using MRI equipment¹.

1 Sherlock, F. G. *Reference Manual for Magnetic Resonance Safety, Implants, and Devices: 2011 Edition*. Biomedical Research Publishing Group, 2011.

PRECAUTIONS: An implant shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Instruments shall be inspected for wear or damage prior to usage. Protect implants against scratching and nicking. Such stress concentrations can lead to failure.

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ADVERSE EFFECTS: Possible adverse effects are pain, discomfort, or abnormal sensations and nerve or soft tissue damage due to the presence of an implant or due to surgical trauma. Fracture of the implant may occur due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion. Implant migration and/or

loosening may occur. Metal sensitivity or histological or allergic or adverse foreign body reaction resulting from implantation of a foreign material may occur. Nerve or soft tissue damage, necrosis of bone or bone resorption, necrosis of the tissue or inadequate healing may result from the presence of an implant or due to surgical trauma.

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Warnings & Precautions

- Any implant contaminated with blood, tissue and/or bodily fluids/matter should be processed according to hospital protocol.
- Do not use an implant if the surface has been damaged. Damaged implants should be discarded.
- Users should wear appropriate personal protective equipment (PPE)
- All users should be qualified personnel with documented evidence of training and competency.

Training should be inclusive of current applicable guidelines and standards and hospital policies.

Manual Processing

Equipment: Soft bristled brush, neutral enzymatic cleaner or neutral detergent with a pH \leq 8.5.

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Ultrasonic Processing

Equipment: Ultrasonic cleaner, neutral enzymatic cleaner or neutral detergent with a pH \leq 8.5. Note: Ultrasonic cleaning may cause additional damage to implants that have surface damage.

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Mechanical Processing

Equipment: Washer, dishwasher, neutral enzymatic cleaner or neutral detergent with a pH \leq 8.5.

CYCLE	MINIMUM TIME (MINUTES)	MINIMUM TEMPERATURE WATER	TYPE OF DETERGENT
Pre-wash	2	Cold tap water	N/A
Enzyme wash	2	Warm tap water	Neutral enzymatic pH \leq 8.5
Wash II	5	Warm tap water (>40°C)	Detergent with pH \leq 8.5
Rinse	2	Warm DI or purified water (>40°C)	N/A
Dry	40	90°C	N/A

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1. Prepare enzymatic and cleaning agents at the use-dilution and temperature recommended by the manufacturer. Fresh solutions should be prepared when existing solutions become grossly contaminated.

2. Place instruments in enzymatic solution until completely submerged. Actuate all moveable parts to allow detergent to contact all surfaces. Soak for a minimum of twenty (20) minutes. Use a nylon soft bristled brush to gently scrub instruments until all visible debris is removed. Pay special attention to hard to reach areas. Pay special attention to any cannulated instruments and clean with an appropriate bottle brush. For exposed springs, coils, or flexible features: Flood the crevices with copious amounts of cleaning solution to flush out any soil. Scrub the surface with a scrub brush to remove all visible soil from the surface and crevices. Bend the flexible area and scrub the surface with a scrub brush. Rotate the part while scrubbing to ensure that all crevices are cleaned.
3. Remove the instruments and rinse thoroughly under running water for a minimum three (3) minutes. Pay special attention to cannulations, and use a syringe to flush any hard to reach areas.
4. Place the instruments, fully submerged, in an ultrasonic unit with cleaning solution. Actuate all moveable parts to allow detergent to contact all

surfaces. Sonicate the instruments for a minimum of ten (10) minutes.

5. Remove the instruments and rinse in deionized water for a minimum of three (3) minutes or until all signs of blood or soil are absent in the rinse stream. Pay special attention to cannulations, and use a syringe to flush any hard to reach areas.
6. Inspect instruments under normal lighting for the removal of visible soil.
7. If visible soil is seen, repeat the sonication and rinse steps above.
8. Remove excess moisture from the instruments with a clean, absorbent, nonshedding wipe.

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1. Prepare enzymatic and cleaning agents at the use-dilution and temperature recommended by the manufacturer. Fresh solutions should be prepared when existing solutions become grossly contaminated.
2. Place instruments in enzymatic solution until completely submerged. Actuate all moveable parts

to allow detergent to contact all surfaces. Soak for a minimum of ten (10) minutes. Use a nylon soft bristled brush to gently scrub instruments until all visible debris is removed. Pay special attention to hard to reach areas. Pay special attention to any cannulated instruments and clean with an appropriate bottle brush.

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- Automated washer/dryer systems are not recommended as the only cleaning method for surgical instruments.

Step	Description
1	Two (2) minute prewash with cold tap water
2	Twenty (20) second enzyme spray with hot tap water
3	One (1) minute enzyme soak
4	Fifteen (15) second cold tap water rinse (x2)
5	Two (2) minute detergent wash with hot tap water (64-66°C/146-150°F)
6	Fifteen (15) second hot tap water rinse
7	Ten (10) second purified water rinse with optional lubricant (64-66°C/146-150°F)
8	Seven (7) minute hot air dry (116°C/240°F)

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- An automated system may be used as a follow up process to manual cleaning.
- Instruments should be thoroughly inspected prior to sterilization to ensure effective cleaning.

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Sterile Product: Sterile product was exposed to a minimum dose of 25.0-kGy gamma irradiation. Resterilization of a device provided sterile should be handled according to hospital protocol. Acumed does not

recommend resterilization of sterile-packaged product.

Non-Sterile Product: Unless clearly labeled as sterile and provided in an unopened sterile package provided by Acumed, all implants and instruments must be considered nonsterile, and sterilized by the hospital prior to use. Nonsterile devices have been validated using the sterilization parameters listed below, in fully loaded trays with all parts placed appropriately.

Sterilization Methods

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

AT-7017, ATF-060, ATM-031, AM-5040, AP-0500, 80-0884, 80-0389, 80-0810, 80-0813, 80-1527/80-1534, 80-0586 These trays have been validated to the parameters below:	
Gravity Displacement Autoclave: NOT RECOMMENDED	
Pre-Vacuum Autoclave:	
Condition:	Wrapped
Exposure Temperature:	270° F (132° C)
Exposure Time:	4 minutes
Dry Time:	30 minutes
AT2-0001, AP-67020, AT2-C700, AT2-M600, AT2-S500 These trays have been validated to the parameters below:	
Gravity Displacement Autoclave: NOT RECOMMENDED	
Pre-Vacuum Autoclave:	
Condition:	Wrapped
Exposure Temperature:	270° F (132° C)
Exposure Time:	10 minutes
Dry Time:	40 minutes

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STORAGE INSTRUCTIONS: Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest first.

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
Practice Settings" og ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

- Lyn-sterilisering anbefales ikke, men hvis den anvendes, bør den kun udføres i følge bestemmelserne i ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

AT-7017, ATF-060, ATM-031, AM-5040, AP-0500, 80-0884, 80-0389, 80-0810, 80-0813, 80-1527/80-1534, 80-0586 Disse bakker er blevet godkendt iht. parametrene angivet herunder:	
Autoklavering, tyngdekraftsforskydning: FRARÅDES	
Autoklavering, prævakuum:	
Tilstand:	Indpakket
Eksponeringstemperatur:	132 °C (270 °F)
Eksponeringsperiode:	4 minutter
Tørreperiode:	30 minutter
AT2-0001, AP-67020, AT2-C700, AT2-M600, AT2-S500 Disse bakker er blevet godkendt iht. parametrene angivet herunder:	
Autoklavering, tyngdekraftsforskydning: FRARÅDES	
Autoklavering, prævakuum:	
Tilstand:	Indpakket
Eksponeringstemperatur:	132 °C (270 °F)
Eksponeringsperiode:	10 minutter
Tørreperiode:	40 minutter

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ADVARSEL: Kun til professionelt brug.

EN ACUTRAK®, ACUTRAK 2® SCREW SYSTEMS

FOR THE PERSONAL ATTENTION OF THE OPERATION SURGEON

DESCRIPTION: Acutrak and Acutrak 2 Screws are designed to provide fixation of various fractures and osteotomies while they heal. The Acutrak and Acutrak 2 Screw Systems include the following implants:

- Acutrak Fusion, Mini, Standard, 4/5, 6/7, and Plus Screws
- AcuTwist Screws
- Acutrak 2 Micro, Mini, Standard, 4.7, 5.5, and 7.5 Screws

The Acutrak and Acutrak 2 Screw Systems also include instruments to facilitate placement of implants.

INFORMATION FOR USE: The surgeon must select the type and size implant that best meets the patient's surgical needs.

Although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be

conveyed to the patient.

SURGICAL TECHNIQUES: Surgical techniques are available describing the uses of this system. It is the responsibility of the surgeon to be familiar with the procedure before use of these products. In addition, it is the responsibility of the surgeon to be familiar with relevant publications and consult with experienced associates regarding the procedure before use.

INDICATIONS: Acutrak and Acutrak 2 screws have the following indications:

Acutrak Fusion, Mini, and Standard, AcuTwist; Acutrak 2 Micro, Mini, Standard, 4.7, and 5.5

Screws: Intended as a fixation device for small bones, bone fragments, and osteotomies. It is not intended for interference or soft tissue fixation.

Acutrak 4/5, 6/7, and Plus; Acutrak 2 7.5 Screws:

This fixation device may be used for fusions, fractures, or osteotomies of the clavicle, humerus, radius, ulna,

ilium, femur, patella, fibula, tibia, talus, malleolus, and calcaneus.

CONTRAINDICATIONS: Contraindications for the system are active or latent infection; sepsis; osteoporosis, insufficient quantity or quality of bone/soft tissue; and material sensitivity. If sensitivity is suspected, tests are performed prior to implantation. Patients who are unwilling or incapable of following postoperative care instructions are contraindicated for these devices. These devices are not intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

WARNINGS: For safe effective use of the implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for the device. The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Device breakage or damage can occur when the implant is subjected to increased loading associated with delayed

union, nonunion, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening, or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant. These cautions include the possibility of the device or treatment failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing, and the possibility of nerve or soft tissue damage related to either surgical trauma or the presence of the implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail. The components of these systems have not been tested for safety, heating, or migration in the MRI environment. The implants may cause distortion and/or block the view of anatomic structures on radiographic images. Similar products have been tested and described in terms of how they may be safely used in post-operative clinical evaluation using MRI equipment.

1 Shellock, F. G. Reference Manual for Magnetic Resonance Safety,

Implants, and Devices: 2011 Edition. Biomedical Research Publishing Group, 2011.

PRECAUTIONS: An implant shall never be reused.

Previous stresses may have created imperfections, which can lead to a device failure. Instruments shall be inspected for wear or damage prior to usage. Protect implants against scratching and nicking. Such stress concentrations can lead to failure.

Mixing implant components from different manufacturers is not recommended for metallurgical, mechanical and functional reasons.

The benefits from implant surgery may not meet the patient's expectations or may deteriorate with time, necessitating revision surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are not uncommon.

ADVERSE EFFECTS: Possible adverse effects are pain, discomfort, or abnormal sensations and nerve or soft tissue damage due to the presence of an implant or due to surgical trauma. Fracture of the implant may occur due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive

force exerted on the implant during insertion.

Implant migration and/or loosening may occur. Metal sensitivity or histological or allergic or adverse foreign body reaction resulting from implantation of a foreign material may occur. Nerve or soft tissue damage, necrosis of bone or bone resorption, necrosis of the tissue or inadequate healing may result from the presence of an implant or due to surgical trauma.

CLEANING:

Implant Cleaning:

Implants should not be re-used. However, any implant that has not been used, but has become soiled, should be processed according to the following:

Warnings & Precautions

- Any implant contaminated with blood, tissue and/or bodily fluids/matter should be processed according to hospital protocol.
- Do not use an implant if the surface has been damaged. Damaged implants should be discarded.
- Users should wear appropriate personal protective equipment (PPE)
- All users should be qualified personnel with

documented evidence of training and competency. Training should be inclusive of current applicable guidelines and standards and hospital policies.

Manual Processing

Equipment: Soft bristled brush, neutral enzymatic cleaner or neutral detergent with a pH ≤ 8.5.

1. Prepare a solution using warm tap water and detergent or cleaner. Follow the enzymatic cleaner or detergent manufacturer's recommendations for use paying close attention to the correct exposure time, temperature, water quality, and concentration.
2. Carefully wash the implant manually. Do not use steel wool or abrasive cleaners on implants.
3. Rinse implant thoroughly with DI or purified water. Use DI or purified water for final rinse.
4. Dry the implant using a clean soft, lint-free cloth to avoid scratching the surface.

Ultrasonic Processing

Equipment: Ultrasonic cleaner, neutral enzymatic cleaner or neutral detergent with a pH ≤ 8.5. Note: Ultrasonic cleaning may cause additional damage to implants that have surface damage.

1. Prepare a solution using warm tap water and

detergent or cleaner. Follow the enzymatic cleaner or detergent manufacturer's recommendations for use paying close attention to the correct exposure time, temperature, water quality, and concentration.

2. Clean implants ultrasonically for a minimum of 15 minutes.
3. Rinse implant thoroughly with DI or purified water. Use DI or purified water for final rinse.
4. Dry the implant using a clean soft, lint-free cloth to avoid scratching the surface.

Mechanical Processing

Equipment: Washer, dishwasher, neutral enzymatic cleaner or neutral detergent with a pH ≤ 8.5.

CYCLE	MINIMUM TIME (MINUTES)	MINIMUM TEMPERATURE WATER	TYPE OF DETERGENT
Pre-wash	2	Cold tap water	N/A
Enzyme wash	2	Warm tap water	Neutral enzymatic pH ≤ 8.5
Wash II	5	Warm tap water (>40°C)	Detergent with pH ≤ 8.5
Rinse	2	Warm DI or purified water (>40°C)	N/A
Dry	40	90°C	N/A

Instrument Cleaning:

Acumed Instruments and Accessories must be thoroughly cleaned before reuse, following the guidelines below.

Warnings & Precautions

- Decontamination of reusable instruments or accessories should occur immediately after completion of the surgical procedure. Do not allow contaminated instruments to dry prior to cleaning/reprocessing. Excess blood or debris should be wiped off to prevent it from drying onto the surface.
- All users should be qualified personnel with documented evidence of training and competency. Training should be inclusive of current applicable guidelines and standards and hospital policies.
- Do not use metal brushes or scouring pads during manual cleaning process.
- Use cleaning agents with low foaming surfactants for manual cleaning in order to see instruments in the cleaning solution. Cleaning agents must be easily rinsed from instruments to prevent residue.
- Mineral oil or silicone lubricants should not be used

on Acumed instruments.

- Neutral pH enzymatic and cleaning agents are recommended for cleaning reusable instruments. It is very important that alkaline cleaning agents are thoroughly neutralized and rinsed from instruments.
- Surgical instruments must be dried thoroughly to prevent rust formation, even if manufactured from high grade stainless steel.
- All instruments must be inspected for cleanliness of surfaces, joints, and lumens, proper function, and wear and tear prior to sterilization.
- Anodized aluminum must not come in contact with certain cleaning or disinfectant solutions. Avoid strong alkaline cleaners and disinfectants or solutions containing iodine, chlorine or certain metal salts. Also, in solutions with pH values above 11, the anodization layer may dissolve.

Manual Cleaning/Disinfection Instructions

1. Prepare enzymatic and cleaning agents at the use-dilution and temperature recommended by the manufacturer. Fresh solutions should be prepared when existing solutions become grossly

contaminated.

2. Place instruments in enzymatic solution until completely submerged. Actuate all moveable parts to allow detergent to contact all surfaces. Soak for a minimum of twenty (20) minutes. Use a nylon soft bristled brush to gently scrub instruments until all visible debris is removed. Pay special attention to hard to reach areas. Pay special attention to any cannulated instruments and clean with an appropriate bottle brush. For exposed springs, coils, or flexible features: Flood the crevices with copious amounts of cleaning solution to flush out any soil. Scrub the surface with a scrub brush to remove all visible soil from the surface and crevices. Bend the flexible area and scrub the surface with a scrub brush. Rotate the part while scrubbing to ensure that all crevices are cleaned.
3. Remove the instruments and rinse thoroughly under running water for a minimum three (3) minutes. Pay special attention to cannulations, and use a syringe to flush any hard to reach areas.

4. Place the instruments, fully submerged, in an ultrasonic unit with cleaning solution. Actuate all moveable parts to allow detergent to contact all surfaces. Sonicate the instruments for a minimum of ten (10) minutes.
5. Remove the instruments and rinse in deionized water for a minimum of three (3) minutes or until all signs of blood or soil are absent in the rinse stream. Pay special attention to cannulations, and use a syringe to flush any hard to reach areas.
6. Inspect instruments under normal lighting for the removal of visible soil.
7. If visible soil is seen, repeat the sonication and rinse steps above.
8. Remove excess moisture from the instruments with a clean, absorbent, nonshedding wipe.

Combination Manual/Automated Cleaning and Disinfecting Instructions

1. Prepare enzymatic and cleaning agents at the use-dilution and temperature recommended by the manufacturer. Fresh solutions should be prepared when existing solutions become grossly

contaminated.

2. Place instruments in enzymatic solution until completely submerged. Actuate all moveable parts to allow detergent to contact all surfaces. Soak for a minimum of ten (10) minutes. Use a nylon soft bristled brush to gently scrub instruments until all visible debris is removed. Pay special attention to hard to reach areas. Pay special attention to any cannulated instruments and clean with an appropriate bottle brush.

Note: Use of a sonicator will aid in thorough cleaning of instruments. Using a syringe or water jet will improve flushing of difficult to reach areas and any closely mated surface.

3. Remove instruments from enzyme solution and rinse in deionized water for a minimum of one (1) minute.
4. Place instruments in a suitable washer/disinfector basket and process through a standard washer/disinfector cycle. The following minimum parameters are essential for thorough cleaning and disinfection.

Step	Description
1	Two (2) minute prewash with cold tap water
2	Twenty (20) second enzyme spray with hot tap water
3	One (1) minute enzyme soak
4	Fifteen (15) second cold tap water rinse (X2)
5	Two (2) minute detergent wash with hot tap water (64-66°C/146-150°F)
6	Fifteen (15) second hot tap water rinse
7	Ten (10) second purified water rinse with optional lubricant (64-66°C/146-150°F)
8	Seven (7) minute hot air dry (116°C/240°F)

Note: Follow washer/disinfector manufacturer's instructions explicitly

Automated Cleaning/Disinfection Instructions

- Automated washer/dryer systems are not recommended as the only cleaning method for surgical instruments.
- An automated system may be used as a follow up process to manual cleaning.
- Instruments should be thoroughly inspected prior to sterilization to ensure effective cleaning.

STERILITY:

System components may be provided sterile or nonsterile.

Sterile Product: Sterile product was exposed to a minimum dose of 25.0-kGy gamma irradiation. Resterilization of a device provided sterile should be handled according to hospital protocol. Acumed does not recommend resterilization of sterile-packaged product.

Non-Sterile Product: Unless clearly labeled as sterile and provided in an unopened sterile package provided by Acumed, all implants and instruments must be considered nonsterile, and sterilized by the hospital prior to use. Nonsterile devices have been validated using the sterilization parameters listed below, in fully-loaded trays with all parts placed appropriately.

Sterilization Methods

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

STORAGE INSTRUCTIONS: Store in a cool dry place and keep away from direct sunlight. Prior to use,

AT-7017, ATF-060, ATM-031, AM-5040, AP-0500, 80-0884, 80-0389, 80-0810, 80-0813, 80-1527/80-1534, 80-0586 These trays have been validated to the parameters below:	
Gravity Displacement Autoclave: NOT RECOMMENDED	
Pre-Vacuum Autoclave:	
Condition:	Wrapped
Exposure Temperature:	270°F (132°C)
Exposure Time:	4 minutes
Dry Time:	30 minutes
AT2-0001, AP-67020, AT2-C000, AT2-M600, AT2-S500 These trays have been validated to the parameters below:	
Gravity Displacement Autoclave: NOT RECOMMENDED	
Pre-Vacuum Autoclave:	
Condition:	Wrapped
Exposure Temperature:	270°F (132°C)
Exposure Time:	10 minutes
Dry Time:	40 minutes

SYMBOL LEGEND	
	Consult instructions for use
	Caution
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Use-by date
	Batch code
	Authorized representative in the European Community
	Manufacturer
	Date of manufacture
	Do not re-sterilize
	Do not re-use
	Upper limit of temperature

OBTAIN THE LATEST VERSION
 Visit www.acumed.net
 for the latest version

inspect product package for signs of tampering, or water contamination. Use oldest first.

APPLICABILITY: These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located.

FURTHER INFORMATION: To request further material, please see the contact information listed on this document.

OBSOLETE
Visit www.acumed.net/ifu
for the latest version.

Cautions: Federal Law (USA) restricts this product sale by or on the order of a physician or hospital. Professional Use Only.

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