ACUMED®

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FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

DESCRIPTION: Acumed Disposable Single Use Sterile Surgical Instruments are available for a wide variety of surgical techniques and implant systems. These instruments should be discarded after each use. The Acumed Bone Bone Graft System instruments are designed to harvest morsellized cancellous bone from the iliac crest, proximal ulna, distal radius, and other donor sites. Acumed Scaphoid Targeting Guide Disks are designed for use in combination with Acumed's Scaphoid Targeting Guide for Acutrak screws in percutaneous fixation of non-displaced fractures of the scaphoid. The disks should be discarded after each use. These instrument systems are disposable, single use instrument and should be discarded after each use.

INFORMATION FOR USE: The instruments are available in various sizes, and the surgeon must select the type and size of implant that best meets the patient's needs and the requirement for the specific orthopedic procedure. The size of the Bone Graft instrument used should meet both the size of the patient and the harvest site. The instruments are not reusable. The type of targeting guide disk is chosen based on the size of implants used.

INDICATIONS: Acumed Single Use Sterile Instruments are designed for use in orthopedic procedures. The Bone Graft System is design for procedures that require autogenous bone graft material. The Scaphoid Targeting Guide Disks are designed for use with Acumed's Scaphoid Targeting Guide in percutaneous Acutrak screw fixation of non-displaced fractures of the Scaphoid bone.

WARNINGS: For safe effective use of any Acumed instrument, the surgeon must be thoroughly familiar with the instrument, the method of application, associated instrumentation, and the recommended surgical technique. Instrument breakage or damage can occur when the instrument is subjected to excessive loads, speeds, or dense bone and care should be taken to avoid such conditions. The patient must be cautioned, preferably in writing, about the use, limitations and possible adverse effects of these type of instruments. The patient must be warned that failure to follow postoperative care instructions can cause the treatment to fail.

PRECAUTIONS: Acumed Disposable Single Use Sterile Surgical Instruments, Acumed Bone Graft System instruments and Acumed Scaphoid Targeting Guide Disks shall never be reused. The instruments should be inspected for defects prior to usage. Previous stresses may have created imperfections, which can lead to a device failure. Protect instruments against scratching and nicking, such stress concentrations can lead to failure. The Acumed Bone Graft System must be used at rpm's that are generated either by hand or with a slow speed drill.

ADVERSE EFFECTS: Use of the Bone Graft System instruments may cause pain, discomfort, or abnormal sensations, nerve damage resulting from the surgical trauma, necrosis, bone resorption, or fracture through the donor site.

STERILITY: This product is provided sterile and was exposed to a minimum dose of 25.0 kGy gamma irradiation. Resterilization may only be performed if the original sterile package has been opened in error using one of the following methods.

Sterilization Methods

Gravity Displacement Autoclave	Full Cycle Time: 270° F (132° C) for 60 minutes Dry Time: 125 minutes
Pre-Vacuum Autoclave	Full Cycle Time: 270° F (132° C) for 25 minutes Dry Time: 70 minutes

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79:2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79:2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

STORAGE INSTRUCTIONS: Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest lots first.

CAUTION: For Professional Use Only.

CAUTION: Federal Law (USA) restricts this product sale by or on the order of a physician or hospital.

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STERILITY: This product is provided sterile and was exposed to a minimum dose of 25.0 kGy gamma irradiation. Resterilization may only be performed if the original sterile package has been opened in error using one of the following methods

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- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

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