ACUMED®

ACU-SINCH REPAIR SYSTEM

MediMark® Europe Sarl.
11 rue Emile ZOLA. BP 2332
38033 GRENOBLE CEDEX 2
FRANCE
+33 4 76 86 43 22

Acumed® LLC
5885 NW Cornelius Pass Road
Hillsboro, OR 97124-9432
+1 (503) 627-9957
www.acumed.net

PKGI-60-B
Effective 6/2012

STERILE EO

REF

LOT

0473

* P K G I - 6 Q - B *

OBSOLETE
Visit www.acumed.net/itd
for the latest version.

ACU-SINCH™ REPARATIONSSYSTEM

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CONTRAINDICATIONS: Contraindications for the system are active or latent infection; sepsis; osteoporosis, insufficient quantity or quality of bone and/or soft tissue; and material sensitivity. If sensitivity is suspected, tests are to be performed prior to implantation. Patients who are unwilling or incapable of following postoperative care instructions are contraindicated for these devices. The device is not indicated as a sole treatment for chronic ligament and tendon injuries. The system may not be suitable for skeletally immature patients, and must not disturb the growth plate.

WARNINGS: For safe effective use of the implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for the device. The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Device breakage (including suture breakage) or damage can occur when the device is subjected to increased loading associated with delayed union, nonunion, or incomplete soft tissue healing. Improper insertion of the device during implantation can increase the possibility of implant loosening and migration, neurovascular injury, or bony fracture. Drilling into or through bone may result in immediate or delayed bone fracture. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant. These cautions include the possibility of the device or treatment failing as a result of loose fixation and/or loosening, stress, excessive activity, or

weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete soft tissue healing; and the possibility of soft tissue damage related to either surgical trauma (i.e. damage to neurovascular structures during preparation such as drilling or implant insertion) or purely due to the presence of an implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail.

MAGNETIC RESONANCE IMAGING (MRI) WARNING: The components of the Acu-Sinch™ Repair System have not been tested for safety, heating, or migration in the MRI environment. The implants may cause distortion and/or block the view of anatomic structures on radiographic images. Similar products have been tested and described in terms of how they may be safely used in post operative clinical evaluation using MRI equipment¹.

¹ Shenoy S. F. Reference Manual for Magnetic Resonance Safety: Implants, and Devices: 2011 Edition. Elsevier Medical Research Publishing Group, 2011.

PRECAUTIONS: An implant shall never be reused. Previous stresses may have created imperfections which can lead to a device failure. Instruments shall be inspected for wear or damage prior to usage. Protect implants against scratching and nicking, as such stress concentrations can lead to failure. Inspect suture for signs of damage or fraying. Prior to use, inspect the product package for signs of tampering, water contamination, and seal integrity.

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- Pain, discomfort or abnormal sensations
- Nerve, vascular or other soft tissue damage
- Necrosis of the tissue or bone
- Bone resorption
- Inadequate healing due to the presence of an implant or due to surgical trauma.

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1 Shellrock, P. en Reference Manual for Magnetic Resonance Safety, Implants, and Devices: 2011 Edition, Biomedical Research Publications, 2011.

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1 Shellcock, F. G. Reference Manual for Magnetic Resonance Safety, Implants, and Devices, 2011 Edition, Biomedical Research Publishing Group, 2011

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ACU-SINCH™ REPAIR SYSTEM

FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

INDICATIONS: The Acu-Sinch™ Repair System is intended to be used in conjunction with a clavicle plate of the Congruent Bone Plate System to provide fixation during the healing of clavicle fractures.

The Acu-Sinch™ Repair System also may be used as a stand-alone system for treatment of acromioclavicular and/or coracoclavicular ligament disruption.

INFORMATION FOR USE: The surgeon must decide based on the injury type and severity what combination of plate and/or soft tissue fixation as well as implantation location is best to meet the patient's surgical needs.

CONTRAINDICATIONS: Contraindications for the system are active or latent infection; sepsis; osteoporosis; insufficient quantity or quality of bone and/or soft tissue; and material sensitivity. If sensitivity is suspected, tests are to be performed prior to implantation. Patients who are unwilling or incapable of following postoperative care instructions are contraindicated for these devices. The device is not indicated as a sole treatment for chronic ligament and tendon injuries. The system may not be suitable for skeletally immature patients, and must not disturb the growth plate.

WARNINGS: For safe effective use of the implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for the device. The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Device breakage (including suture breakage) or damage can occur when the device is subjected to increased loading associated with delayed union, nonunion, or incomplete soft tissue healing. Improper insertion of the device during implantation can increase the possibility of implant loosening and migration, neurovascular injury, or bony fracture. Drilling into or through bone may result in immediate or delayed bone fracture. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant. These cautions include the possibility of the device or treatment failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete soft tissue healing; and the

possibly of soft tissue damage related to either surgical trauma (i.e. damage to neurovascular structures during preparation such as drilling or implant insertion) or purely due to the presence of an implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail.

MAGNETIC RESONANCE IMAGING (MRI) WARNING: The components of the Acu-Sinch™ Repair System have not been tested for safety, heating, or migration in the MRI environment. The implants may cause distortion and/or block the view of anatomic structures on radiographic images. Similar products have been tested and described in terms of how they may be safely used in post-operative clinical evaluation using MRI equipment¹.

¹ Shellock, F.G. Reference Manual to Magnetic Resonance Safety, Implants, and Devices, 2011 Edition. Biomedical Research Publishing Group, 2011.

PRECAUTIONS: An implant shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Instruments shall be inspected for wear or damage prior to usage. Protect implants against scratching and nicking, as such stress concentrations can lead to failure. Inspect suture for signs of damage or fraying. Prior to use, inspect the product package for signs of tampering, water contamination, and seal integrity.

ADVERSE EFFECTS: Possible adverse affects are:

- Pain, discomfort, or abnormal sensations
- Nerve, vascular, or other soft tissue damage
- Necrosis of the tissue or bone
- Bone resorption
- Inadequate healing due to the presence of an implant or due to surgical trauma.

Implant migration and/or loosening, and suture breakage may occur due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the device during insertion and during the healing period, as well as inadequate postoperative protocol or failed patient compliance to the postoperative protocol. Implant material sensitivity or histological or allergic reaction resulting from implantation of a foreign material may occur.

IMPLANT MATERIAL SPECIFICATIONS:

Implants: The metal implants are made out of implantable grade titanium.
Suture: The suture is composed of Ultra High Molecular Weight Polyethylene (UHMWPE).

STERILITY: Implants and instruments are provided sterile. Sterility was achieved by using the ethylene oxide gas sterilization method to Sterility Assurance Level (SAL) of 10⁻⁶. Do not re-sterilize.

STORAGE INSTRUCTIONS: Store in a cool dry place and keep away from direct sunlight. Do not expose suture to heat. Use oldest first.

CAUTION: Federal Law (USA) restricts this product sale by or on the order of a physician or hospital.