# ACUMED®

## ACUMED® CANNULATED SCREW SYSTEM

MediMark® Europe Sarl.

11 rue Emile ZOLA. BP 2332  
38033 GRENOBLE CEDEX 2  
FRANCE  
+33 4 76 86 43 22

EC REP

Acumed®

5885 NW Cornelius Pass Road  
Hillsboro, OR 97124-9432  
+1 (503) 627-9957  
www.acumed.net

PKG-66-A  
EFFECTIVE 6/2013

STERILE R

REF

LOT

0473

\* P K G I - 6 6 - A \*

Visit [www.acumed.net/ru](http://www.acumed.net/ru) for the latest version.

## TABLE OF CONTENTS

Click on language to navigate to page.

English – US	3
Dansk – DA	8
English – EN	14
Deutsch – DE	19
Ελληνικά – EL	25
Español – ES	31
Français – FR	37
Italiano – IT	43
Nederlands – NL	49
Norsk – NO	55
Português – PT	60
Suomi – FI	66
Svenska – SV	72
Türkçe – TR	77

**OBSOLETE**  
Visit [www.acumed.net/ifu](http://www.acumed.net/ifu)  
for the latest version.

## US ACUMED® CANNULATED SCREW SYSTEM

### FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

**DESCRIPTION:** The Acumed Cannulated Screws are designed to provide fixation of various fractures, fusions, and osteotomies while they heal.

The Acumed Cannulated Screw System also includes instruments to facilitate placement of implants.

**INFORMATION FOR USE:** The surgeon must select the type and size implant that best meets the patient's surgical needs.

Although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be conveyed to the patient.

**SURGICAL TECHNIQUES:** Surgical techniques are available describing the uses of this system. It is the responsibility of the surgeon to be familiar with the procedure before use of these products. In addition, it is the responsibility of the surgeon to be familiar with

relevant publications and consult with experienced associates regarding the procedure before use.

**INDICATIONS:** The Acumed Cannulated Screw System is generally intended for fixation of fractures, fusion and osteotomies of large and small bones appropriate for size of device, which may include the following:

Minimally invasive reconstruction of fractures and joints; Adjuvant for osteosynthesis in complex joint fractures; Multifragment joint fractures; Simple metaphyseal fractures; Fractures of the wrist, ankle, elbow, and shoulder; Condylar fractures; Epiphyseal and metaphyseal fractures in children; Osteochondritis dissecans; Osteo-Chondral Fractures; Ligament avulsion injuries; Ligament fixation; Other small fragment, cancellous bone fractures; Small joint fusion; Areas where accurate screw placement is vital; Metatarsal and phalangeal osteotomies; Fractures of the tarsals, metatarsals and other fractures of the foot; Avulsion fractures and fractures of metatarsal V; Tarso-metatarsal and metatarso-phalangeal arthrodesis;

Tarsal fusions; Calcaneal and talar fractures; Subtalar arthrodesis; Ankle arthrodesis; Fractures of small joints, such as ankle fractures and navicular fractures; Fractures of the fibula, malleolus, and calcaneus; Distal tibia and pilon fractures; Acetabular fractures; Other fractures of the pelvic ring; Sacroiliac joint disruptions; Fractures of the femoral head and neck; Supracondylar femoral fractures; Slipped capital femoral epiphyses; An adjunct to DHS in basilar neck fractures; Pediatric femoral neck fractures; Intercondylar femur fractures; Intracapsular fractures of the hip; Fractures of the distal femur and proximal tibia; Patellar fractures; Tibial plateau fractures; Small fragments of the hand and wrist; Fractures of the carpals and metacarpals; Carpal and metacarpal arthrodesis; Scaphoid fracture and other fractures of the hand; Phalangeal and interphalangeal fractures; Fractures of the ulna and radius; Radial head fractures; Fractures of the olecranon and distal humerus; Humeral head fractures; Ligament fixation at the proximal humerus;

and Glenoid fractures.

Washers may be used with the screws in certain applications.

**CONTRAINDICATIONS:** Contraindications for the system are active or latent infection, sepsis, osteoporosis, insufficient quantity or quality of bone/ soft tissue, and material sensitivity. If sensitivity is suspected, tests are performed prior to implantation. Patients who are unwilling or incapable of following postoperative care instructions are contraindicated for these devices. These devices are not intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

**WARNINGS:** For safe effective use of this implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for this device. The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Device breakage or damage can occur when the implant is subjected to increased loading associated with delayed union, nonunion, or incomplete healing.

Improper insertion of the device during implantation can increase the possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant including the possibility of the device or treatment failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing including possible nerve or soft tissue damage related to either surgical trauma or the presence of the implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail. The implant system has not been evaluated for safety and compatibility in the MR (Magnetic Resonance) environment, nor has it been tested for heating or migration in the MR environment.

**PRECAUTIONS:** An implant shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Instruments shall be inspected for wear or damage prior to usage. Protect implants against scratching and nicking. Such stress

concentrations can lead to failure.

Mixing implant components from different manufacturers is not recommended for metallurgical, mechanical and functional reasons.

The benefits from implant surgery may not meet the patient's expectations or may deteriorate with time, necessitating revision surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are not uncommon.

**ADVERSE EFFECTS:** Possible adverse effects are pain, discomfort, or abnormal sensations and nerve or soft tissue damage due to the presence of an implant or due to surgical trauma. Fracture of the implant may occur due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion. Implant migration and/ or loosening may occur. Metal sensitivity or histological or allergic reaction resulting from implantation of a foreign material may occur. Nerve or soft tissue damage, necrosis of bone or bone resorption, necrosis of the tissue or inadequate healing may result from the presence of an implant or due to surgical trauma.

## CLEANING:

**Implant Cleaning:** Implants should not be re-used. However, any implant that has not been used, but has become soiled, should be processed according to the following:

### Warnings & Precautions

- Any implant contaminated with blood, tissue and/or bodily fluids/matter should be processed according to hospital protocol.
- Do not use an implant if the surface has been damaged. Damaged implants should be discarded.
- Users should wear appropriate personal protective equipment (PPE)
- All users should be qualified personnel with documented evidence of training and competency. Training should be inclusive of current applicable guidelines and standards and hospital policies.

### Manual Processing

*Equipment: Soft bristled brush, neutral enzymatic cleaner or neutral detergent with a pH < 8.5.*

1. Prepare a solution using warm tap water and

detergent or cleaner. Follow the enzymatic cleaner or detergent manufacturer's recommendations for use paying close attention to the correct exposure time, temperature, water quality, and concentration.

2. Carefully wash the implant manually. Do not use steel wool or abrasive cleaners on implants.
3. Rinse implant thoroughly with DI or purified water. Use DI or purified water for final rinse.
4. Dry the implant using a clean soft, lint-free cloth to avoid scratching the surface.

### Ultrasonic Processing

*Equipment: Ultrasonic cleaner, neutral enzymatic cleaner or neutral detergent with a pH < 8.5. Note: Ultrasonic cleaning may cause additional damage to implants that have surface damage.*

1. Prepare a solution using warm tap water and detergent or cleaner. Follow the enzymatic cleaner or detergent manufacturer's recommendations for use paying close attention to the correct exposure time, temperature, water quality, and concentration.

2. Clean implants ultrasonically for a minimum of 15 minutes.
3. Rinse implant thoroughly with DI or purified water. Use DI or purified water for final rinse.
4. Dry the implant using a clean soft, lint-free cloth to avoid scratching the surface.

### Mechanical Processing

*Equipment: Washer/disinfecter, neutral enzymatic cleaner or neutral detergent with a pH < 8.5.*

CYCLE	MINIMUM TIME (MINUTES)	MINIMUM TEMPERATURE/WATER	TYPE OF DETERGENT
Pre-wash	2	Cold tap water	N/A
Enzyme wash	2	Warm tap water	Neutral enzymatic pH ≤ 8.5
Wash II	5	Warm tap water (>40°C)	Detergent with pH ≤ 8.5
Rinse	2	Warm DI or purified water (>40°C)	N/A
Dry	40	90°C	N/A

**Instrument Cleaning Requirements:** Acumed Instruments and Accessories must be thoroughly cleaned before reuse. Decontamination of reusable instruments or accessories should occur immediately after completion of the surgical procedure. Excess blood or debris should be wiped off to prevent it from drying onto the surface. Using an enzymatic cleaning product such as Enzol, prepare the cleaning solution according to manufacturer's instructions. Avoiding prolonged soaking, submerge the device in the enzymatic cleaning solution for 5 minutes. After 5 minutes, use a soft bristle brush to remove any detergent residue. Instruments with lumens also should be held under water (larger side up) ensuring the lumen is flushed with water. Place the device in an ultrasonic cleaner filled with deionized water, ensuring that there is no air trapped in any of the crevices by gently agitating the device. Sonicate the device for 5 minutes. Remove the device from the ultrasonic cleaner and rinse in water by submerging it completely. Gently wipe the device with a clean, dry cloth, then air dry.

Even surgical instruments manufactured from high

grade stainless steel must be dried thoroughly to prevent rust formation, and all devices must be inspected for cleanliness of surfaces, joints, and lumens, proper function, and wear and tear prior to sterilization.

Instruments manufactured from different metals should be processed separately to avoid electrolytic action between the different metals.

#### *Exposed Springs, Coils, or Flexible Features*

Flood the crevices with copious amounts of cleaning solution to flush out any soil. Scrub the surface with a scrub brush to remove all visible soil from the surface and crevices. Bend the flexible area and scrub the surface with a scrub brush. Rotate the part while scrubbing to ensure that all crevices are cleaned.

#### *Anodized Aluminum*

Anodized aluminum must not come in contact with certain cleaning or disinfectant solutions. Avoid strong alkaline cleaners and disinfectants or solutions containing iodine, chlorine or certain metal salts. Also, in solutions with pH values above 11, the anodization layer may dissolve.

#### **STERILITY:**

System components may be provided sterile or nonsterile.

**Sterile Product:** Sterile product was exposed to a minimum dose of 25.0-kGy gamma irradiation. Resterilization of a device provided sterile should be handled according to hospital protocol. Acumed does not recommend resterilization of sterile-packaged product.

**Non-Sterile Product:** Unless clearly labeled as sterile and provided in an unopened sterile package provided by Acumed, all implants and instruments must be considered nonsterile, and sterilized by the hospital prior to use. Nonsterile devices have been validated using the sterilization parameters listed below, in fully-loaded trays with all parts placed appropriately.

#### *Sterilization Methods*

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and

ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

80-1823/80-1824: This tray has been validated to the parameters below:	
Gravity Displacement Autoclave: NOT RECOMMENDED	
Pre-Vacuum Autoclave:	
Condition:	Wrapped
Exposure Temperature:	270° F (132° C)
Exposure Time:	4 minutes
Dry Time:	30 minutes

**STORAGE INSTRUCTIONS:** Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest lots first.

SYMBOL LEGEND	
	Consult instructions for use
	Caution
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Use-by date
	Batch code
	Authorized representative in the European Community
	Manufacturer
	Date of manufacture
	Do not re-sterilize
	Do not re-use
	Upper limit of temperature

**APPLICABILITY:** These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located.

**FURTHER INFORMATION:** To request further material, please see the contact information listed on this document.

Cautions: Federal Law (USA) restricts this product sale by or on the order of a physician or hospital. Professional Use Only.

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Minimally invasive reconstruction of fractures and joints; Adjuvant for osteosynthesis in complex joint fractures; Multifragment joint fractures; Simple metaphyseal fractures; Fractures of the wrist, ankle, elbow, and shoulder; Condylar fractures; Epiphyseal and metaphyseal fractures in children; Osteochondritis dissecans, Osteo-Chondral Fractures; Ligament avulsion injuries; Ligament fixation; Other small fragment, cancellous bone fractures; Small joint fusion; Areas where accurate screw placement is vital; Metatarsal and phalangeal osteotomies; Fractures of the tarsals, metatarsals and other fractures of the foot; Avulsion fractures and fractures of metatarsal V; Tarso-

metatarsal and metatarso-phalangeal arthrodesis; Tarsal fusions; Calcaneal and talar fractures; Subtalar arthrodesis; Ankle arthrodesis; Fractures of small joints, such as Ankle fractures and Navicular fractures; Fractures of the fibula, malleolus, and calcaneus; Distal tibia and pilon fractures; Acetabular fractures; Other fractures of the pelvic ring; Sacroiliac joint disruptions; Fractures of the femoral head and neck; Supracondylar femoral fractures; Slipped capital femoral epiphyses; An adjunct to DHS in basilar neck fractures; Pediatric femoral neck fractures; Intercondylar femur fractures; Intracapsular fractures of the hip; Fractures of the distal femur and proximal tibia; Patellar fractures; Tibial plateau fractures; Small fragments of the hand and wrist; Fractures of the carpals and metacarpals; Carpal and metacarpal arthrodesis; Scaphoid fracture and other fractures of the hand; Phalangeal and interphalangeal fractures; Fractures of the ulna and radius; Radial head fractures; Fractures of the olecranon and distal humerus; Humeral head

fractures; Ligament fixation at the proximal humerus; and Glenoid fractures.

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**WARNINGS:** For safe effective use of this implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for this device. The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Device breakage or damage can occur when the implant is subjected to increased loading associated

with delayed union, nonunion, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant including the possibility of the device or treatment failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing including possible nerve or soft tissue damage related to either surgical trauma or the presence of the implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail. The implant system has not been evaluated for safety and compatibility in the MR (Magnetic Resonance) environment, nor has it been tested for heating or migration in the MR environment.

**PRECAUTIONS:** An implant shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Instruments shall be inspected for wear or damage prior to usage. Protect

implants against scratching and nicking. Such stress concentrations can lead to failure.

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**ADVERSE EFFECTS:** Possible adverse effects are pain, discomfort, or abnormal sensations and nerve or soft tissue damage due to the presence of an implant or due to surgical trauma. Fracture of the implant may occur due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion. Implant migration and/or loosening may occur. Metal sensitivity or histological or allergic reaction resulting from implantation of a foreign material may occur. Nerve or soft tissue damage, necrosis of bone or bone resorption, necrosis of the tissue or inadequate healing

may result from the presence of an implant or due to surgical trauma.

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### Manual Processing

*Equipment: Soft bristled brush, neutral enzymatic cleaner or neutral detergent with a pH ≤ 8.5.*

1. Prepare a solution using warm tap water and detergent or cleaner. Follow the enzymatic cleaner or detergent manufacturer's recommendations for use paying close attention to the correct exposure time, temperature, water quality, and concentration.
2. Carefully wash the implant manually. Do not use steel wool or abrasive cleaners on implants.
3. Rinse implant thoroughly with DI or purified water. Use DI or purified water for final rinse.
4. Dry the implant using a clean soft, lint-free cloth to avoid scratching the surface.

### Ultrasonic Processing

*Equipment: Ultrasonic cleaner, neutral enzymatic cleaner or neutral detergent with a pH ≤ 8.5. Note: Ultrasonic cleaning may cause additional damage to implants that have surface damage.*

1. Prepare a solution using warm tap water and detergent or cleaner. Follow the enzymatic cleaner or detergent manufacturer's recommendations for use paying close attention to the correct exposure time, temperature, water quality, and

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### Mechanical Processing

*Equipment: Washer/disinfector, neutral enzymatic cleaner or neutral detergent with a pH ≤ 8.5.*

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**Instrument Cleaning Requirements:** Acumed Instruments and Accessories must be thoroughly cleaned before reuse. Decontamination of reusable

instruments or accessories should occur immediately after completion of the surgical procedure. Excess blood or debris should be wiped off to prevent it from drying onto the surface. Using an enzymatic cleaning product such as Enzol, prepare the cleaning solution according to manufacturer's instructions. Avoiding prolonged soaking, submerge the device in the enzymatic cleaning solution for 5 minutes. After 5 minutes, use a soft bristle brush to remove any detergent residue. Instruments with lumens also should be held under water (larger side up) ensuring the lumen is flushed with water. Place the device in an ultrasonic cleaner filled with deionized water, ensuring that there is no air trapped in any of the crevices by gently agitating the device. Sonicate the device for 5 minutes. Remove the device from the ultrasonic cleaner and rinse in water by submerging it completely. Gently wipe the device with a clean, dry cloth, then air dry.

Even surgical instruments manufactured from high grade stainless steel must be dried thoroughly to prevent rust formation, and all devices must be inspected for cleanliness of surfaces, joints, and

lumens, proper function, and wear and tear prior to sterilization.

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#### *Sterilization Methods*

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

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**APPLICABILITY:** These materials contain information about products that may or may not be available in any particular country or may be available under

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different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located.

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