 **acumed**<sup>®</sup>  
**Biotrak™ Resorbable Fixation System**

**OBSOLETE**  
Visit [www.acumed.net/ifu](http://www.acumed.net/ifu)  
for the latest version.

MediMark® Europe Sarl.  
11 rue Emile ZOLA. BP 2332  
38033 GRENOBLE CEDEX 2  
FRANCE  
+33 4 76 86 43 22

Acumed® LLC  
5885 NE Cornelius Pass Road  
Hillsboro, OR 97124-9432  
+1.503.627.9957  
acumed.net

**77013001-M**  
**EFFECTIVE 03-2018**

## TABLE OF CONTENTS

Click on a language to navigate to the page

English – US .....	3
Dansk – DA .....	7
English – EN .....	11
Deutsch – DE .....	15
Ελληνικά – EL .....	19
Español – ES .....	23
Français – FR .....	27
Italiano – IT .....	31
Nederlands – NL.....	35
Norsk – NO .....	39
Portugues – PT .....	43
Suomi – FI .....	47
Svenska – SV .....	51
Türkçe – TR .....	55

OBSELETE  
Visit [www.acumed.net/ifu](http://www.acumed.net/ifu)  
for the latest version.

## BIOTRAK™ RESORBABLE FIXATION SYSTEM

### FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

**DESCRIPTION:** The ACUMED resorbable fixation devices are designed to provide fixation of specific fractures, osteotomies, and arthrodeses while they heal.

**INDICATIONS:** This fixation device is indicated for small bone fractures, osteotomies, and arthrodeses, cancellous fragments, and osteochondral fragments in the upper and lower extremities.

**CONTRAINDICATIONS:** This device is contraindicated in the presence of active or latent infection, sepsis, osteoporosis, insufficient quantity and/or quality of bone. This device is not intended for use in any indication not covered above.

**IMPLANT MATERIAL SPECIFICATIONS:** The implants are made of PLLA (Poly-L-lactic acid).

**SURGICAL INSTRUMENT MATERIAL SPECIFICATIONS:** The instruments are made of various grades of stainless steel, titanium, aluminum, and polymers evaluated for biocompatibility.

**IMPLANT INFORMATION FOR USE:** The surgeon must select the type and size of implant that best meets the patient's surgical needs. Although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be conveyed to the patient.

#### SURGICAL INSTRUMENT INFORMATION FOR USE:

Instruments provided with this system may be single use or reusable.

- The user must refer to the instrument's label to determine whether the instrument is single use or reusable. Single use instruments are labeled with a "do not re-use" symbol as described in the Symbol Legend section, below.
- Single use instruments must be discarded after a single use.
- Reusable instruments have a limited lifespan. Prior to and after each use, reusable instruments must be inspected where applicable for sharpness, wear, damage, proper cleaning, corrosion and integrity of the connecting mechanisms. Particular care should be paid to drivers, drill bits and instruments used for cutting or implant insertion.

**SURGICAL TECHNIQUES:** Surgical techniques are available describing the uses of this system. It is the responsibility of the surgeon to be familiar with the procedure before use of these

products. In addition, it is the responsibility of the surgeon to be familiar with relevant publications and consult with experienced associates regarding the procedure before use. Surgical Techniques can be found on the Acumed website (acumed.net).

**IMPLANT WARNINGS:** For safe and effective use of this implant, the surgeon must be thoroughly familiar with the implant, material, methods of application, instruments, and recommended surgical technique for this device. The device is not designed to withstand the stress of weight bearing, load bearing, and/or excessive activity. Device breakage and/or damage can occur when the implant is subjected to increased loading associated with delayed union, non-union, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant including the possibility of the device failing due to aforementioned causes. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail. In addition, because of the thermal sensitivity of bioabsorbable materials, the device shall not be utilized if the dot in the middle of the temperature sticker has turned black. A black dot on the sensor signifies that the environmental temperature may have exceeded the softening temperature of the bioabsorbable material during storage and/or transit. Exceeding the softening temperature of the material can lead to degradation of the mechanical properties and/or warping.

**SURGICAL INSTRUMENT WARNINGS:** For safe effective use of any Acumed instrument, the surgeon must be familiar with the instrument, the method of application, and the recommended surgical technique. Instrument breakage or damage, as well as tissue damage, can occur when an instrument is subjected to excessive loads, excessive speeds, dense bone, improper use or unintended use. The patient must be cautioned, preferably in writing as to the risks associated with these types of instruments.

**IMPLANT PRECAUTIONS:** An implant shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Protect implants against scratching and nicking. Such stress concentrations can lead to failure. Bending plates multiple times may weaken the device and could lead to premature implant fracture and failure. Mixing implant components from different manufacturers is not recommended for metallurgical, mechanical and functional reasons. The benefits from implant surgery may not meet the patient's expectations or may deteriorate with time, necessitating revision surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are not uncommon.

**SURGICAL INSTRUMENT PRECAUTIONS:** Single use surgical instruments shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Protect instruments against scratching and nicking, such stress concentrations can lead to failure.

**ADVERSE EFFECTS:** Fracture of the implant due to excessive loading, incomplete or inadequate healing, or excessive force during insertion; Implant migration and/or loosening; Sterile inflammation as a result of a body reaction to the degradation products of the absorbable material; Pain, discomfort, or abnormal sensations due to the presence of an implant; Nerve damage resulting from surgical trauma; Bone necrosis or bone resorption.

the center has not turned black. If the temperature of the packaging has exceeded the softening temperature of the material during its lifetime, the central dot on the temperature sensor will be black and this device must not be implanted. Use the oldest lots first because the bioabsorbable material has a finite shelf life.

#### CLEANING:

This product is provided sterile, and should not be re-cleaned.

#### STERILITY:

This product is only provided sterile.

**Sterile Product:** Sterile product was exposed to a minimum dose of 25.0-kGy gamma irradiation. Acumed does not recommend re-sterilization of sterile-packaged product. If sterile packaging is damaged, the incident must be reported to Acumed. The product must not be used, and must be returned to Acumed.

**STORAGE INSTRUCTIONS:** Store implants in a cool dry place and keep away from direct sunlight. The implant should be stored in a temperature and humidity controlled environment below 120° F. Prior to use, inspect packaging for signs of tampering or water contamination. Check the expiration date on the box to ensure that the shelf life of the product has not been exceeded. Also, please inspect the temperature sensor on the box to ensure that

**APPLICABILITY:** These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located.

**FURTHER INFORMATION:** To request further material, please see the contact information listed on this document.

SYMBOL LEGEND	
	Consult instructions for use
	Caution
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Use-by date
	Catalogue number
	Batch code
	Authorized representative in the European Community
	Manufacturer
	Date of manufacture
	Do not re-sterilize
	Do not re-use
	Upper limit of temperature

**Cautions: Federal Law (USA) restricts this product sale by or on the order of a physician or hospital. Professional Use Only.**

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#### **SURGICAL INSTRUMENT INFORMATION FOR USE:**

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**IMPLANT PRECAUTIONS:** An implant shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Protect implants against scratching and nicking. Such stress concentrations can lead to failure. The benefits from implant surgery may not meet the patient's expectations or may deteriorate with time, necessitating revision surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are not uncommon.

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**ADVERSE EFFECTS:** Fracture of the implant due to excessive loading, incomplete or inadequate healing, or excessive force during insertion; Implant migration and/or loosening; Sterile inflammation as a result of a body reaction to the degradation products of the absorbable material; Pain, discomfort, or abnormal

sensations due to the presence of an implant; Nerve damage resulting from surgical trauma; Bone necrosis or bone resorption.

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**Sterile Product:** Sterile product was exposed to a minimum dose of 25.0-kGy gamma irradiation. Acumed does not recommend resterilization of sterile-packaged product. If sterile packaging is damaged, the incident must be reported to Acumed. The product must not be used, and must be returned to Acumed.

**STORAGE INSTRUCTIONS:** Store implants in a cool dry place and keep away from direct sunlight. The implant should be stored in a temperature and humidity controlled environment below 120° F. Prior to use, inspect packaging for signs of tampering or water contamination. Check the expiration date on the box to ensure that the shelf life of the product has not been exceeded. Also, please inspect the temperature sensor on the box to ensure that the center has not turned black. If the temperature of the packaging has exceeded the softening temperature of the material during its lifetime, the central dot on the temperature sensor will be black and this device must not be implanted. Use the oldest lots first because the bioabsorbable material has a finite shelf life.

**STORAGE INSTRUCTIONS:** Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest lots first.

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