ACUMED®

ACUMED BONE GRAFT SYSTEM

MediMark® Europe Sarl.
11 rue Emile ZOLA. BP 2332
38033 GRENOBLE CEDEX 2
FRANCE
+33 4 76 86 43 22

Acumed®
5885 NW Cornelius Pass Road
Hillsboro, OR 97124-9432
+1 (503) 627-9957
www.acumed.net

PKGI-23-B
Effective 11/2009

0473

* P K G I - 2 3 - B *

OBSOLETE
Visit www.acumed.net/itd
for the latest version

ACUMEDS KNOGLEGRAFTSSYSTEM

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ACUMED BONE GRAFT SYSTEM

FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

DESCRIPTION: Acumed bone graft instruments are designed to harvest morselized cancellous bone from the iliac crest, distal radius, and distal femur.

INFORMATION FOR USE: The instruments are available in various sizes, and it is important to choose the instrument that meets both the size of the patient and the harvesting site. The instruments are reusable; however, they have a limited life span. Prior to each use, the instruments must be inspected for sharpness and integrity of the connecting mechanisms.

INDICATIONS: Procedures which require autogenous cancellous bone graft material.

CONTRAINDICATIONS: Active or latent infection. Osteoporosis, insufficient quantity or quality of bone/soft tissue. Material sensitivity. If suspected, tests are performed prior to implantation. Sepsis. Patients who are unwilling or incapable of following postoperative care instructions.

WARNINGS: For safe and effective use of this instrument, the surgeon must be thoroughly familiar with the instrument, the method of application, and the recommended surgical technique. Instrument breakage or damage can occur when the instrument is subjected to excessive loads, speeds, or extremely dense bone. Improper insertion of the device can lead to soft tissue and/or bony complications. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this instrument. The patient must be warned that failure to follow postoperative care instructions can cause the treatment to fail.

PRECAUTIONS: The instrument should be inspected for wear or damage prior to each usage. The instrument must be used at rpm's which are generated either by hand or with a slow speed drill.

ADVERSE EFFECTS: Pain, discomfort, or abnormal sensations due to the presence of an implant. Nerve damage resulting from the surgical trauma. Necrosis of bone or bone resorption. Necrosis of the tissue or inadequate healing.

STERILITY: This product is provided non-sterile. Sterilization may be performed using one of the following methods below.

Sterilization Methods

Gravity Displacement Autoclave	270° F (132° C) for 30 minutes
Pre-Vacuum Autoclave	270° F (132° C) for 4 minutes

- Please consider your equipment manufacturer's written instructions for the specific sterilizer and load configuration being used and current AORN standards and recommended practices.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

STORAGE INSTRUCTIONS: Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest lots first.

CAUTION: For Professional Use Only.

CAUTION: Federal Law (USA) restricts this product sale by or on the order of a physician or hospital.

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DESCRIPTION: The Acumed Bone Graft instruments are designed to harvest morselized cancellous bone from the iliac crest, distal radius, and distant femur.

INFORMATION FOR USE: The instruments are available in various sizes, and it is important to choose the instrument that meets both the size of the patient and the harvest site. The instruments are reusable; however, they have a limited life span. Prior to each use, the instruments must be inspected for sharpness and integrity of connecting mechanisms.

INDICATIONS: Procedures that require autogenous bone graft material.

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WARNINGS: For safe effective use of this implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for this device. Instrument breakage or damage can occur when the instrument is subjected to excessive loads, speeds or extremely dense bone. Improper insertion of the device can lead to soft tissue and/or bony complications. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this instrument. The patient must be warned that failure to follow postoperative care instructions can cause the treatment to fail.

PRECAUTIONS: The instruments should be inspected for wear or damage prior to usage. The instrument must be used at rpm's that are generated either by hand or with a slow speed drill.

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- Please consider your equipment manufacturer's written instructions for the specific sterilizer and load configuration being used and current AORN standards and recommended practices.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST19: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

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