 **acumed**[®]
Single Use Instruments

CE
0086

MediMark[®] Europe Sarl.
11 rue Emile ZOLA, BP 2332
38033 GRENOBLE CEDEX 2
FRANCE
+33 4 76 86 43 22

EC REP

Acumed[®] LLC
5885 NW Cornelius Pass Road
Hillsboro, OR 97124-9432
+1.503.627.9957
acumed.net

OBSOLETE
Visit www.acumed.net/ifu
for the latest version.

PKGI-57-F
EFFECTIVE 08-2016

TABLE OF CONTENTS

Click on a language to navigate to the page

| | |
|-----------------------|----|
| English – US | 3 |
| Dansk – DA | 9 |
| English – EN | 15 |
| Deutsch – DE | 21 |
| Ελληνικά – EL | 27 |
| Español – ES | 33 |
| Français – FR | 39 |
| Italiano – IT | 45 |
| Nederlands – NL | 51 |
| Norsk – NO | 57 |
| Português – PT | 63 |
| Suomi – FI | 69 |
| Svenska – SV | 75 |
| Türkçe – TR | 81 |

OBSOLETE
Visit www.acumed.net/ifu
for the latest version.

ACUMED® SINGLE USE INSTRUMENTS

FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

DESCRIPTION: Acumed Disposable, Single Use Surgical Instruments are available for a wide variety of surgical techniques and implant systems. Acumed Disposable, Single Use Surgical Tap Instruments are only available for use with Acumed bone plates. These instruments should be discarded after each use.

SURGICAL INSTRUMENT INFORMATION FOR USE:

Instruments provided with this system may be single use or reusable.

- The user must refer to the instrument's label to determine whether the instrument is single use or reusable. Single use instruments are labeled with a "do not re-use" symbol as described in the Symbol Legend section, below.
- Single use instruments must be discarded after a single use.
- Reusable instruments have a limited lifespan. Prior to and after each use, reusable instruments must be inspected where applicable for sharpness, wear, damage, proper cleaning, corrosion and integrity of the connecting mechanisms. Particular care should be paid to drivers, drill bits and instruments used for cutting or implant insertion.

SURGICAL INSTRUMENT WARNINGS: For safe effective use of any Acumed instrument, the surgeon must be familiar with the instrument, the method of application, and the recommended surgical technique. Instrument breakage or damage, as well as tissue damage, can occur when an instrument is subjected to excessive loads, excessive speeds, dense bone, improper use or unintended use. The patient must be cautioned, preferably in writing as to the risks associated with these types of instruments. The surgical tap should only be used with Acumed bone plates as recommended in the surgical technique. Use with plates or implants not recommended in the surgical technique can cause plate, tap, and/or treatment failure. The patient must be cautioned, preferably in writing, about the risks associated with these types of instruments.

SURGICAL INSTRUMENT PRECAUTIONS: Single use surgical instruments shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Protect instruments against scratching and nicking, such stress concentrations can lead to failure.

Tapping a plate using a plate tap will cause titanium debris to be generated, which should be removed. Failure to remove the plate debris can cause, among other complications, inflammation, cartilage damage, and patient discomfort. The taps are single surgery use and should be discarded after each

surgery or if the tap becomes dull or damaged. If the resistance increases while using a tap, discard the tap immediately. Breakage to the tap can occur due to excessive torque or levering and care should be taken to avoid such conditions. Should breakage occur, carefully remove all tap pieces.

CLEANING:

Implant Cleaning: Cleaning instructions for a particular implant can be reviewed in the Instructions for Use assigned to the implant part number.

Instrument Cleaning: Acumed Instruments and Accessories must be thoroughly cleaned before reuse, following the guidelines below:

Warnings & Precautions

- Decontamination of reusable instruments or accessories should occur immediately after completion of the surgical procedure. Do not allow contaminated instruments to dry prior to cleaning/ reprocessing. Excess blood or debris should be wiped off to prevent it from drying onto the surface.
- All users should be qualified personnel with documented evidence of training and competency. Training should be inclusive of current applicable guidelines and standards and hospital policies.
- Do not use metal brushes or scouring pads during manual cleaning process.
- Use cleaning agents with low foaming surfactants for manual cleaning in order to see instruments in the cleaning solution.

Cleaning agents must be easily rinsed from instruments to prevent residue.

- Mineral oil or silicone lubricants should not be used on Acumed instruments.
- Neutral pH enzymatic and cleaning agents are recommended for cleaning reusable instruments. It is very important that alkaline cleaning agents are thoroughly neutralized and rinsed from instruments.
- Surgical instruments must be dried thoroughly to prevent rust formation, even if manufactured from high grade stainless steel.
- All instruments must be inspected for cleanliness of surfaces, joints, and lumens, proper function, and wear and tear prior to sterilization.
- Anodized aluminum must not come in contact with certain cleaning or disinfectant solutions. Avoid strong alkaline cleaners and disinfectants or solutions containing iodine, chlorine or certain metal salts. Also, in solutions with pH values above 11, the anodization layer may dissolve.
- To prolong the life of torque limiting drivers, moving parts must be lubricated with a commercial water-based surgical grade instrument lubricant (instrument milk) after every use. Lubrication should be performed after cleaning and prior to sterilization. Follow lubricant manufacturer's instructions.

Manual Cleaning/Disinfection Instructions

1. Prepare enzymatic and cleaning agents at the use-dilution and temperature recommended by the manufacturer. Fresh solutions should be prepared when existing solutions become grossly contaminated.

2. Place instruments in enzymatic solution until completely submerged. Actuate all moveable parts to allow detergent to contact all surfaces. Soak for a minimum of twenty (20) minutes. Use a nylon soft bristled brush to gently scrub instruments until all visible debris is removed. Pay special attention to hard to reach areas. Pay special attention to any cannulated instruments and clean with an appropriate bottle brush. For exposed springs, coils, or flexible features: Flood the crevices with copious amounts of cleaning solution to flush out any soil. Scrub the surface with a scrub brush to remove all visible soil from the surface and crevices. Bend the flexible area and scrub the surface with a scrub brush. Rotate the part while scrubbing to ensure that all crevices are cleaned.
3. Remove the instruments and rinse thoroughly under running water for a minimum three (3) minutes. Pay special attention to cannulations, and use a syringe to flush any hard to reach areas.
4. Place the instruments, fully submerged, in an ultrasonic unit with cleaning solution. Actuate all moveable parts to allow detergent to contact all surfaces. Sonicate the instruments for a minimum of ten (10) minutes.
5. Remove the instruments and rinse in deionized water for a minimum of three (3) minutes or until all signs of blood or soil are absent in the rinse stream. Pay special attention to cannulations, and use a syringe to flush any hard to reach areas.
6. Inspect instruments under normal lighting for the removal of visible soil.

7. If visible soil is seen, repeat the sonication and rinse steps above.
8. Remove excess moisture from the instruments with a clean, absorbent, nonshedding wipe.

Combination Manual/Automated Cleaning and Disinfecting Instructions

1. Prepare enzymatic and cleaning agents at the use-dilution and temperature recommended by the manufacturer. Fresh solutions should be prepared when existing solutions become grossly contaminated.
2. Place instruments in enzymatic solution until completely submerged. Actuate all moveable parts to allow detergent to contact all surfaces. Soak for a minimum of ten (10) minutes. Use a nylon soft bristled brush to gently scrub instruments until all visible debris is removed. Pay special attention to hard to reach areas. Pay special attention to any cannulated instruments and clean with an appropriate bottle brush. *Note: Use of a sonicator will aid in thorough cleaning of instruments. Using a syringe or water jet will improve flushing of difficult to reach areas and any closely mated surface.*
3. Remove instruments from enzyme solution and rinse in deionized water for a minimum of one (1) minute.
4. Place instruments in a suitable washer/ disinfecter basket and process through a standard washer/disinfecter cycle. The following minimum parameters are essential for thorough cleaning and disinfection.

| Step | Description |
|------|--|
| 1 | Two (2) minute prewash with cold tap water |
| 2 | Twenty (20) second enzyme spray with hot tap water |
| 3 | One (1) minute enzyme soak |
| 4 | Fifteen (15) second cold tap water rinse (X2) |
| 5 | Two (2) minute detergent wash with hot tap water (64–66°C/146–150°F) |
| 6 | Fifteen (15) second hot tap water rinse |
| 7 | Ten (10) second purified water rinse with optional lubricant (64–66°C/146–150°F) |
| 8 | Seven (7) minute hot air dry (116°C/240°F) |

Note: Follow washer/disinfectant manufacturer's instructions explicitly

Automated Cleaning/Disinfection Instructions

- Automated washer/dryer systems are not recommended as the only cleaning method for surgical instruments.
- An automated system may be used as a follow up process to manual cleaning.
- Instruments should be thoroughly inspected prior to sterilization to ensure effective cleaning.

STERILITY:

System components may be provided sterile or nonsterile.

Sterile Product: Sterile product was exposed to a minimum dose of 25.0-kGy gamma irradiation. Acumed does not recommend resterilization of sterile-packaged product. If sterile packaging is damaged, the incident must be reported to Acumed. The product must not be used, and must be returned to Acumed.

Non-Sterile Product: Unless clearly labeled as sterile and provided in an unopened sterile package provided by Acumed, all implants and instruments must be considered nonsterile, and sterilized by the hospital prior to use. Sterilization of nonsterile devices has been validated using the sterilization parameters listed below, where devices are provided in fully-loaded trays with all parts placed appropriately.

Sterilization Methods

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79:

2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

| | |
|---|-----------------|
| Gravity Displacement Autoclave: NOT RECOMMENDED | |
| Pre-Vacuum Autoclave: | |
| Condition: | Wrapped |
| Exposure Temperature: | 270° F (132° C) |
| Exposure Time: | 4 minutes |
| Dry Time: | 30 minutes |

OBSOLETE
Visit www.acumed.net/ifu
for the latest version.

STORAGE INSTRUCTIONS: Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest lots first.

APPLICABILITY: These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located.

FURTHER INFORMATION: To request further material, please see the contact information listed on this document.

| SYMBOL LEGEND | |
|--|----------------------------------|
|  | Consult instructions for use |
|  | Caution |
|  | Sterilized using ethylene oxide |
|  | Sterilized using irradiation |
|  | Use-by date |
|  | Catalogue number |
|  | Batch code |
|  | Authorized representative in the |
|  | Manufacturer |
|  | Date of manufacture |
|  | Do not re-sterilize |
|  | Do not re-use |
|  | Upper limit of temperature |

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician. For Professional Use Only.

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Tapping a plate using a plate tap will cause titanium debris to be generated, which should be removed. Failure to remove the plate debris can cause, among other complications, inflammation, cartilage damage, and patient discomfort. The

taps are single surgery use and should be discarded after each surgery or if the tap becomes dull or damaged. If the resistance increases while using a tap, discard the tap immediately. Breakage to the tap can occur due to excessive torque or levering and care should be taken to avoid such conditions. Should breakage occur, carefully remove all tap pieces.

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Warnings & Precautions

- Decontamination of reusable instruments or accessories should occur immediately after completion of the surgical procedure. Do not allow contaminated instruments to dry prior to cleaning/ reprocessing. Excess blood or debris should be wiped off to prevent it from drying onto the surface.
- All users should be qualified personnel with documented evidence of training and competency. Training should be inclusive of current applicable guidelines and standards and hospital policies.
- Do not use metal brushes or scouring pads during manual cleaning process.

- Use cleaning agents with low foaming surfactants for manual cleaning in order to see instruments in the cleaning solution. Cleaning agents must be easily rinsed from instruments to prevent residue.
- Mineral oil or silicone lubricants should not be used on Acumed instruments.
- Neutral pH enzymatic and cleaning agents are recommended for cleaning reusable instruments. It is very important that alkaline cleaning agents are thoroughly neutralized and rinsed from instruments.
- Surgical instruments must be dried thoroughly to prevent rust formation, even if manufactured from high grade stainless steel.
- All instruments must be inspected for cleanliness of surfaces, joints, and lumens, proper function, and wear and tear prior to sterilization.
- Anodized aluminum must not come in contact with certain cleaning or disinfectant solutions. Avoid strong alkaline cleaners and disinfectants or solutions containing iodine, chlorine or certain metal salts. Also, in solutions with pH values above 11, the anodization layer may dissolve.
- To prolong the life of torque limiting drivers, moving parts must be lubricated with a commercial water-based surgical grade instrument lubricant (instrument milk) after every use. Lubrication should be performed after cleaning and prior to sterilization. Follow lubricant manufacturer's instructions.

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1. Prepare enzymatic and cleaning agents at the use-dilution and temperature recommended by the manufacturer. Fresh

- solutions should be prepared when existing solutions become grossly contaminated.
- Place instruments in enzymatic solution until completely submerged. Actuate all moveable parts to allow detergent to contact all surfaces. Soak for a minimum of twenty (20) minutes. Use a nylon soft bristled brush to gently scrub instruments until all visible debris is removed. Pay special attention to hard to reach areas. Pay special attention to any cannulated instruments and clean with an appropriate bottle brush. For exposed springs, coils, or flexible features: Flood the crevices with copious amounts of cleaning solution to flush out any soil. Scrub the surface with a scrub brush to remove all visible soil from the surface and crevices. Bend the flexible area and scrub the surface with a scrub brush. Rotate the part while scrubbing to ensure that all crevices are cleaned.
 - Remove the instruments and rinse thoroughly under running water for a minimum three (3) minutes. Pay special attention to cannulations, and use a syringe to flush any hard to reach areas.
 - Place the instruments, fully submerged, in an ultrasonic unit with cleaning solution. Actuate all moveable parts to allow detergent to contact all surfaces. Sonicate the instruments for a minimum of ten (10) minutes.
 - Remove the instruments and rinse in deionized water for a minimum of three (3) minutes or until all signs of blood or soil are absent in the rinse stream. Pay special attention to cannulations, and use a syringe to flush any hard to reach areas.
 - Inspect instruments under normal lighting for the removal of visible soil.
 - If visible soil is seen, repeat the sonication and rinse steps above.
 - Remove excess moisture from the instruments with a clean, absorbent, nonshedding wipe.
- Combination Manual/Automated Cleaning and Disinfecting Instructions**
- Prepare enzymatic and cleaning agents at the use-dilution and temperature recommended by the manufacturer. Fresh solutions should be prepared when existing solutions become grossly contaminated.
 - Place instruments in enzymatic solution until completely submerged. Actuate all moveable parts to allow detergent to contact all surfaces. Soak for a minimum of ten (10) minutes. Use a nylon soft bristled brush to gently scrub instruments until all visible debris is removed. Pay special attention to hard to reach areas. Pay special attention to any cannulated instruments and clean with an appropriate bottle brush. *Note: Use of a sonicator will aid in thorough cleaning of instruments. Using a syringe or water jet will improve flushing of difficult to reach areas and any closely mated surface.*
 - Remove instruments from enzyme solution and rinse in deionized water for a minimum of one (1) minute.
 - Place instruments in a suitable washer/ disinfecter basket and process through a standard washer/disinfecter cycle. The following minimum parameters are essential for thorough cleaning and disinfection.

| Step | Description |
|--|--|
| 1 | Two (2) minute prewash with cold tap water |
| 2 | Twenty (20) second enzyme spray with hot tap water |
| 3 | One (1) minute enzyme soak |
| 4 | Fifteen (15) second cold tap water rinse (X2) |
| 5 | Two (2) minute detergent wash with hot tap water (64–66°C/146–150°F) |
| 6 | Fifteen (15) second hot tap water rinse |
| 7 | Ten (10) second purified water rinse with optional lubricant |
| 8 | Seven (7) minute hot air dry (116°C/240°F) |
| <i>Note: Follow washer/disinfectant manufacturer's instructions explicitly</i> | |

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| SYMBOL LEGEND | |
|---|----------------------------------|
|  | Consult instructions for use |
|  | Caution |
|  | Sterilized using ethylene oxide |
|  | Sterilized using irradiation |
|  | Use-by date |
|  | Catalogue number |
|  | Batch code |
|  | Authorized representative in the |
|  | Manufacturer |
|  | Date of manufacture |
|  | Do not re-sterilize |
|  | Do not re-use |
|  | Upper limit of temperature |

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