ACUMED®

ACUMED® CALC-JAK™ CALCANEAL INSTRUMENTATION SYSTEM

OBSOLETE
Visit www.acumed.net/ifu
for the latest version.

MediMark® Europe Sarl.
11 rue Emile ZOLA. BP 2332
38033 GRENOBLE CEDEX 2
FRANCE
+33 4 76 86 43 22

EC REP

Acumed®
5885 NW Cornelius Pass Road
Hillsboro, OR 97124-9432
+1 (503) 627-9957

www.acumed.net

PKG P65-A
EFFECTIVE 7/2013

REF LOT ! CE

* P K G I - 6 5 - A *

TABLE OF CONTENTS

Click on language to navigate to page.

English – US	3
Dansk – DA	7
English – EN	11
Deutsch – DE	15
Ελληνικά – EL	20
Español – ES	25
Français – FR	30
Italiano – IT	34
Nederlands – NL	38
Norsk – NO	43
Português – PT	47
Suomi – FI	52
Svenska – SV	56
Türkçe – TR	60

OBSOLETE
Visit www.acumed.net/ifu
for the latest version.

US ACUMED® CALC-JAK™ CALCANEAL INSTRUMENTATION SYSTEM

FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

DESCRIPTION: The CALC-JAK Calcaneal Instrumentation System is designed to assist implantation of orthopedic implants indicated for calcaneal fractures.

INFORMATION FOR USE: The surgeon should select the type and size of implant that best meets the patient's requirements, according to the Instructions for Use assigned to the implant component part number.

Although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be conveyed to the patient.

SURGICAL TECHNIQUES: Surgical techniques are available describing use of this instrument in conjunction with Acumed systems. It is the responsibility of the surgeon to be familiar with the

procedure before use of these products. In addition, it is the responsibility of the surgeon to be familiar with relevant publications and consult with experienced associates regarding the procedure before use.

INDICATIONS: The CALC-JAK Calcaneal Instrumentation System is designed to assist implantation of orthopedic implants indicated for calcaneal fractures.

Specific implant indications are given in the Instructions for Use assigned to the implant component part number, and should be reviewed prior to use.

CONTRAINDICATIONS: Specific implant contraindications are given in the Instructions for Use assigned to the implant component part number, and should be reviewed prior to use.

WARNINGS: For safe effective use of any Acumed instrument, the surgeon must be familiar with the instrument, the method of application, and the

recommended surgical technique. The patient must be cautioned, preferably in writing, as to the risks associated with these types of instruments.

PRECAUTIONS: Instrument breakage or damage can occur when an instrument is subjected to excessive loads, speeds, or dense bone. The use of an instrument for tasks other than those for which it is indicated may result in damaged or broken instruments. Proper care and inspection after each use is critical; while the instrumentation provided is reusable, it does have a limited life span. Prior to and after each use, the instruments must be inspected where applicable for sharpness, wear, damage, proper cleaning, corrosion and integrity of the connecting mechanisms. Particular care should be paid to hex drivers, drill bits and instruments used for cutting or implant insertion.

ADVERSE EFFECTS: Specific implant adverse effects are given in the Instructions for Use assigned to the implant component part number, and should be reviewed prior to use.

CLEANING

Implant Cleaning Requirements: This system does not contain implants. Cleaning instructions for a particular implant can be reviewed in the Instructions for Use assigned to the implant part number.

Instrument Cleaning Requirements: Acumed Instruments and Accessories must be thoroughly cleaned before reuse, following the guidelines below.

Warnings & Precautions

- Decontamination of reusable instruments or accessories should occur immediately after completion of the surgical procedure. Do not allow contaminated instruments to dry prior to cleaning/reprocessing. Excess blood or debris should be wiped off to prevent it from drying onto the surface.
- All users should be qualified personnel with documented evidence of training and competency. Training should be inclusive of current applicable guidelines and standards and hospital policies.
- Do not use metal brushes or scouring pads during manual cleaning process.
- Use cleaning agents with low foaming surfactants for manual cleaning in order to see instruments in the cleaning solution. Cleaning agents must be easily rinsed from instruments to prevent residue.
- Mineral oil or silicone lubricants should not be used on Acumed instruments.
- Neutral pH enzymatic and cleaning agents are recommended for cleaning reusable instruments. It is very important that alkaline cleaning agents are thoroughly neutralized and rinsed from instruments.
- Surgical instruments must be dried thoroughly to prevent rust formation, even if manufactured from high grade stainless steel.
- All instruments must be inspected for cleanliness of surfaces, joints, and lumens, proper function, and wear and tear prior to sterilization.
- Anodized aluminum must not come in contact with certain cleaning or disinfectant solutions. Avoid strong alkaline cleaners and disinfectants or solutions containing iodine, chlorine or certain metal salts. Also, in solutions with pH values above 11, the anodization layer may dissolve.

Manual Cleaning/Disinfection Instructions

1. Prepare enzymatic and cleaning agents at the use-dilution and temperature recommended by the manufacturer. Fresh solutions should be prepared when existing solutions become grossly contaminated.
2. Place instruments in enzymatic solution until completely submerged. Actuate all moveable parts to allow detergent to contact all surfaces. Soak for a minimum of twenty (20) minutes. Use a nylon soft bristled brush to gently scrub instruments until all visible debris is removed. Pay special attention to hard to reach areas. Pay special attention to any cannulated instruments and clean with an appropriate bottle brush. For exposed springs, coils, or flexible features: Flood the crevices with copious amounts of cleaning solution to flush out any soil. Scrub the surface with a scrub brush to remove all visible soil from the surface and crevices. Bend the flexible area and scrub the surface with a scrub brush. Rotate the part while scrubbing to ensure that all crevices are cleaned.
3. Remove the instruments and rinse thoroughly

under running water for a minimum three (3) minutes. Pay special attention to cannulations, and use a syringe to flush any hard to reach areas.

- Place the instruments, fully submerged, in an ultrasonic unit with cleaning solution. Actuate all moveable parts to allow detergent to contact all surfaces. Sonicate the instruments for a minimum of ten (10) minutes.
- Remove the instruments and rinse in deionized water for a minimum of three (3) minutes or until all signs of blood or soil are absent in the rinse stream. Pay special attention to cannulations, and use a syringe to flush any hard to reach areas.
- Inspect instruments under normal lighting for the removal of visible soil.
- If visible soil is seen, repeat the sonication and rinse steps above.
- Remove excess moisture from the instruments with a clean, absorbent, nonshedding wipe.

Automated Cleaning/Disinfection Instructions

- Automated washer/dryer systems are not recommended as the only cleaning method for

surgical instruments.

- An automated system may be used as a follow up process to manual cleaning.
- Instruments should be thoroughly inspected prior to sterilization to ensure effective cleaning.

STERILITY:

System components are provided nonsterile. All instruments must be sterilized by the hospital prior to use. Instruments have been validated using the sterilization parameters listed below, in fully-loaded trays with all parts placed appropriately.

Sterilization Methods

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2010 — Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements

of ANSI/AAMI ST79: 2010 — Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

80-1408/80-1409: Calcaneal Instrumentation Tray Lid/Base	
These trays have been validated to the parameters below:	
Gravity Displacement Autoclave: NOT RECOMMENDED	
Pre-Vacuum Autoclave:	
Condition:	Wrapped
Exposure Temperature:	270° F (132° C)
Exposure Time:	4 minutes
Dry Time:	30 minutes

STORAGE INSTRUCTIONS: Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest lots first.

SYMBOL LEGEND	
	Consult instructions for use
	Caution
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Use-by date
	Batch code
	Authorized representative in the European Community
	Manufacturer
	Date of manufacture
	Do not re-sterilize
	Do not re-use
	Upper limit of temperature

APPLICABILITY: These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries.

Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located.

FURTHER INFORMATION: To request further material, please see the contact information listed on this document.

OBSOLETE
 Visit www.aclumed.net/ifu for the latest version.

Cautions: Federal Law (USA) restricts this product sale by or on the order of a physician or hospital. Professional Use Only.

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Warnings & Precautions

- Decontamination of reusable instruments or accessories should occur immediately after completion of the surgical procedure. Do not allow contaminated instruments to dry prior to cleaning/reprocessing. Excess blood or debris should be wiped off to prevent it from drying onto the surface.
- All users should be qualified personnel with documented evidence of training and competency. Training should be inclusive of current applicable guidelines and standards and hospital policies.
- Do not use metal brushes or scouring pads during manual cleaning process.
- Use cleaning agents with low foaming surfactants

for manual cleaning in order to see instruments in the cleaning solution. Cleaning agents must be easily rinsed from instruments to prevent residue.

- Mineral oil or silicone lubricants should not be used on Acumed instruments.
- Neutral pH enzymatic and cleaning agents are recommended for cleaning reusable instruments. It is very important that alkaline cleaning agents are thoroughly neutralized and rinsed from instruments.
- Surgical instruments must be dried thoroughly to prevent rust formation, even if manufactured from high grade stainless steel.
- All instruments must be inspected for cleanliness of surfaces, joints, and lumens, proper function, and wear and tear prior to sterilization.
- Anodized aluminum must not come in contact with certain cleaning or disinfectant solutions. Avoid strong alkaline cleaners and disinfectants or solutions containing iodine, chlorine or certain metal salts. Also, in solutions with pH values above 11, the anodization layer may dissolve.

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4. Place the instruments, fully submerged, in an ultrasonic unit with cleaning solution. Actuate all moveable parts to allow detergent to contact all surfaces. Sonicate the instruments for a minimum of ten (10) minutes.

5. Remove the instruments and rinse in deionized water for a minimum of three (3) minutes or until all signs of blood or soil are absent in the rinse stream. Pay special attention to cannulations, and use a syringe to flush any hard to reach areas.

6. Inspect instruments under normal lighting for the removal of visible soil.

7. If visible soil is seen, repeat the sonication and rinse steps above.

8. Remove excess moisture from the instruments with a clean, absorbent, nonshedding wipe.

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Sterilization Methods

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements

of ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

80-1823/80-1824: This tray has been validated to the parameters below:	
Gravity Displacement Autoclave: NOT RECOMMENDED	
Pre-Vacuum Autoclave:	
Condition:	Wrapped
Exposure Temperature:	270° F (132° C)
Exposure Time:	4 minutes
Dry Time:	30 minutes

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OBSOLETE
 Visit www.aelumed.net/ifu for the latest version.

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