## Acu-Loc® Wrist Spanning Plate

MediMark® Europe Sarl.  
11 rue Emile ZOLA. BP 2332  
38033 GRENOBLE CEDEX 2  
FRANCE  
+33 4 76 86 43 22

EC REP

Acumed®  
5885 NW Cornelius Pass Road  
Hillsboro, OR 97124-9432  
+1 (503) 627-9957  
[www.acumed.net](http://www.acumed.net)

PKGI-70-C  
EFFECTIVE 6/2014

STERILE R

REF

LOT

0473

\* P K G I - 7 Q - C \*

OBSOLETE  
Visit [www.acumed.net/ifu](http://www.acumed.net/ifu)  
for the latest version.

# TABLE OF CONTENTS

Click on language to navigate to page.

English – US .....	3
Dansk – DA .....	6
English – EN .....	9
Deutsch – DE .....	12
Ελληνικά – EL .....	15
Español – ES .....	18
Français – FR .....	21
Italiano – IT .....	24
Nederlands – NL .....	27
Norsk – NO .....	30
Português – PT .....	33
Suomi – FI .....	36
Svenska – SV .....	39
Türkçe – TR .....	42

**OBSOLETE**  
Visit [www.acumed.net/ifu](http://www.acumed.net/ifu)  
for the latest version.

# US ACUMED® ACU-LOC® WRIST SPANNING PLATE

## FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

**DESCRIPTION:** The Acu-Loc Wrist Spanning Plate is designed to provide temporary fixation for the wrist. The Acu-Loc Wrist Spanning Plate uses the instruments provided with the Acu-Loc 2 system to facilitate placement of implants.

**INFORMATION FOR USE:** Physiological dimensions limit the sizes of implant appliances. The surgeon must select the type and size that best meets the patient's requirements for close adaptation and firm seating with adequate support.

Although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be conveyed to the patient.

**SURGICAL TECHNIQUES:** Surgical techniques are available describing the uses of this system. It is the responsibility of the surgeon to be familiar with the procedure before use of these products.

In addition, it is the responsibility of the surgeon to be familiar with relevant publications and consult with experienced associates regarding the procedure before use.

**INDICATIONS:** The Acu-Loc Wrist Spanning Plate is indicated for fixation of fractures, osteotomies, and non-unions of the radius.

**CONTRAINDICATIONS:** Contraindications for the system are active or latent infection, sepsis, osteoporosis, insufficient quantity or quality of bone or soft tissue and material sensitivity. If material sensitivity is suspected, tests should be performed prior to implantation. Patients who are unwilling or incapable of following postoperative care instructions are contraindicated for these devices. These devices are not intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

**WARNINGS:** For safe effective use of the implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the

recommended surgical technique for the device. The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Device breakage or damage can occur when the implant is subjected to increased loading associated with delayed union, nonunion, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant. These cautions include the possibility of the device or treatment failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing, and the possibility of nerve or soft tissue damage related to either surgical trauma or the presence of the implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/

or treatment to fail. The components of these systems have not been tested for safety, heating, or migration in the MRI environment. The implants may cause distortion and/or block the view of anatomic structures on radiographic images. Similar products have been tested and described in terms of how they may be safely used in post-operative clinical evaluation using MRI equipment<sup>1</sup>.

*1 Shellock, F. G. Reference Manual for Magnetic Resonance Safety, Implants, and Devices: 2011 Edition. Biomedical Research Publishing Group, 2011.*

**PRECAUTIONS:** An implant shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Instruments shall be inspected for wear or damage prior to usage. Protect implants against scratching and nicking. Such stress concentrations can lead to failure.

Mixing implant components from different manufacturers is not recommended for metallurgical, mechanical and functional reasons.

The benefits from implant surgery may not meet the patient's expectations or may deteriorate with time, necessitating revision

surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are not uncommon.

**ADVERSE EFFECTS:** Possible adverse effects are pain, discomfort, or abnormal sensations and nerve or soft tissue damage due to the presence of an implant or due to surgical trauma. Fracture of the implant may occur due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion. Implant migration and/or loosening may occur. Metal sensitivity, histological, allergic or adverse foreign body reaction resulting from implantation of a foreign material may occur. Nerve or soft tissue damage, necrosis of bone or bone resorption, necrosis of the tissue or inadequate healing may result from the presence of an implant or due to surgical trauma.

**CLEANING:** This product is provided sterile, and should not be re-cleaned.

**STERILITY:** This product is only provided sterile. Sterile product was exposed to a minimum dose of 25.0-kGy gamma irradiation. Resterilization of a device

provided sterile should be handled according to hospital protocol. Acumed does not recommend resterilization of sterile-packaged product.

**STORAGE INSTRUCTIONS:** Store in a cool, dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest lots first.

SYMBOL LEGEND	
	Consult instructions for use
	Caution
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Use-by date
	Catalogue number
	Batch code
	Authorized representative in the European Community
	Manufacturer
	Date of manufacture
	Do not re-sterilize
	Do not re-use
	Upper limit of temperature

**APPLICABILITY:** These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located.

**FURTHER INFORMATION:** To request further material, please see the contact information listed on this document.

Visit [www.acumed.net/ifu](http://www.acumed.net/ifu) for the latest version.

**Cautions: Federal Law (USA) restricts this product sale by or on the order of a physician or hospital. Professional Use Only.**

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