 **acumed**<sup>®</sup>  
Suture Passer

MediMark<sup>®</sup> Europe Sarl.  
11 rue Emile ZOLA. BP 2332  
38033 GRENOBLE CEDEX 2  
FRANCE  
+33 4 76 86 43 22

EC REP

Acumed<sup>®</sup>  
5885 NW Cornelius Pass Road  
Hillsboro, OR 97124-9432  
+1 (503) 627-9957  
[www.acumed.net](http://www.acumed.net)

PKGI-74-A  
EFFECTIVE 2/2014

STERILE EO

REF

LOT

CE

0473

\* P K G I - 7 4 - A \*

OBSOLETE  
Visit [www.acumed.net/ifu](http://www.acumed.net/ifu)  
for the latest version.

# TABLE OF CONTENTS

Click on language to navigate to page.

English – US .....	3
Dansk – DA .....	5
English – EN .....	7
Deutsch – DE .....	9
Ελληνικά – EL .....	12
Español – ES .....	14
Français – FR .....	16
Italiano – IT .....	18
Nederlands – NL .....	20
Norsk – NO .....	23
Português – PT .....	25
Suomi – FI .....	27
Svenska – SV .....	29
Türkçe – TR .....	31

**OBSOLETE**  
Visit [www.acumed.net/ifu](http://www.acumed.net/ifu)  
for the latest version.

# US ACUMED® SUTURE PASSER

## FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

**DESCRIPTION:** The Suture Passer is a sterile single-use surgical instrument.

### **MATERIAL SPECIFICATIONS:**

*Handle:* ABS 348.

*Cannula:* 304 Stainless Steel.

*Loop:* Nylon.

### **INFORMATION FOR USE:**

1. Prior to use, inspect the product package for seal integrity and signs of tampering or water contamination.
2. The Suture Passer shall be inspected for damage prior to use.
3. Insert the Suture Passer through a bony socket or tunnel or implant hole or tissue under direct visualization.
4. Feed the surgical suture through the Suture Passer loop.
5. Pull the suture with the Suture Passer through a bony socket or tunnel or implant hole or tissue under direct visualization.
6. Discard the Suture Passer after the

surgery. The Suture Passer is not reusable.

Although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be conveyed to the patient.

**SURGICAL TECHNIQUES:** It is the responsibility of the surgeon to be familiar with the procedure before use of these products. In addition, it is the responsibility of the surgeon to be familiar with relevant publications and consult with experienced associates regarding the procedure before use. Surgical Techniques can be found on the Acumed website ([acumed.net](http://acumed.net))

**INDICATIONS:** The Suture Passer instrument is indicated for passing surgical suture through bony sockets, tunnels, implant holes and tissue in orthopedic surgery.

**CONTRAINDICATIONS:** The Suture Passer instrument is not indicated for implantation.

**WARNINGS:** The Suture Passer must be

used under direct visualization. Extreme care should be used during insertion to avoid unintentional puncture of any internal organs. The surgeon must be thoroughly familiar with the instrument, the method of application, associated instrumentation, and the recommended surgical technique. Instrument breakage or damage can occur when the instrument is subjected to excessive loads, excessively bent, or used in holes that are too small for the Suture Passer to easily pass suture strands through and care should be taken to avoid such conditions.

**PRECAUTIONS:** The Suture Passer is intended for single use only. Re-sterilization or reuse of this device may result in compromised device performance including device failure and other unintended damage.

**ADVERSE EFFECTS:** Possible adverse effects are pain, discomfort, or abnormal sensations and nerve or soft tissue damage due to the presence of an implant or due to surgical trauma.

**CLEANING:** This product is provided sterile, and should not be re-cleaned.

**STERILITY:** This product is only provided sterile. Sterility was achieved by using the ethylene oxide gas sterilization method to Sterility Assurance Level (SAL) of  $10^{-6}$ . Do not re-sterilize.

**STORAGE INSTRUCTIONS:** Store in a cool dry place and keep away from direct sunlight. Do not expose to heat. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest lots first.

SYMBOL LEGEND	
	Consult instructions for use
	Caution
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Use-by date
	Catalogue number
	Batch code
	Authorized representative in the European Community
	Manufacturer
	Date of manufacture
	Do not re-sterilize
	Do not re-use
	Upper limit of temperature

**APPLICABILITY:** These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located.

**FURTHER INFORMATION:** To request further material, please see the contact information listed on this document.

**Cautions: Federal Law (USA) restricts this product sale by or on the order of a physician or hospital. Professional Use Only.**

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**ADVARSEL: Kun til professionelt brug.**

FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

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**SURGICAL TECHNIQUES:** It is the responsibility of the surgeon to be familiar with the procedure before use of these products. In addition, it is the responsibility of the surgeon to be familiar with relevant publications and consult with experienced associates regarding the procedure before use. Surgical techniques can be found on the Acumed website ([acumed.net](http://acumed.net))

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**WARNINGS:** The Suture Passer must be used under direct visualization. Extreme care should be used during insertion to avoid unintentional puncture of any internal organs. The surgeon must be thoroughly familiar with the instrument, the method of application, associated instrumentation, and the recommended surgical technique. Instrument breakage or damage can occur when the instrument is subjected to excessive loads, excessively bent, or used in holes that are too small for the Suture Passer to easily pass suture strands through and care should be taken to avoid such conditions.

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**ADVERSE EFFECTS:** Possible adverse effects are pain, discomfort, or abnormal sensations and nerve or soft tissue damage due to the presence of an implant

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*Cannula:* 304 Stainless Steel.

*Loop:* Nylon.

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1. Prior to use, inspect the product package for seal integrity and signs of tampering or water contamination.
2. The Suture Passer shall be inspected for damage prior to use.
3. Insert the Suture Passer through a bony socket or tunnel or implant hole or tissue under direct visualization.
4. Feed the surgical suture through the Suture Passer loop.
5. Pull the suture with the Suture Passer through a bony socket or tunnel or implant hole or tissue under direct visualization.

6. Discard of the Suture Passer after the surgery. The Suture Passer is not reusable.

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