



Bimalleolar Ankle Fracture Dislocation Treatment Using Plate and Screw Fixation

Case Study

Kent Ellington, MD

A 57-year-old female with a closed bimalleolar fracture dislocation was treated with the Acumed Ankle Plating System 3 and Small Fragment Base Set.



Pre-op and post-op x-rays

Acumed® is a global leader of innovative orthopaedic and medical solutions.



We are dedicated to developing products, service methods, and approaches that improve patient care.

Case Study | Kent Ellington, MD



Patient History

Patient is a 57-year-old who tripped while walking her dog in August 2015. She sustained a closed bimalleolar fracture dislocation that was reduced in the emergency room and splinted. She was presented to the clinic two days later and went to the operating room on post-injury day six.

Treatment

Treatment included open reduction internal fixation (ORIF) of her ankle fracture. The Acumed Ankle Plating System 3 and Small Fragment Base Set were used to provide fixation. This required a routine lateral and medial approach to the ankle. The patient was placed into a supine position. The oblique fibula was reduced and fixed initially with a lag screw then with a neutralization one-third tubular plate. Next, the medial malleolar fragment was reduced then fixed with two 4.0 mm cannulated partially threaded screws with excellent compression.

Discussion

This is a common fracture pattern. The necessary tools, instruments, plates, and screws required to complete ORIF of ankle fractures are included in the Acumed Ankle Plating System 3. This patient has regained full range of motion and had no pain at final follow-up.



6 months postoperative



LEX70-07-A | Effective: 2017/04 | © 2017 Acumed® LLC

Acumed® Headquarters
5885 NW Cornelius Pass Road
Hillsboro, OR 97124
Office: 888.627.9957
Fax: 503.520.9618
www.acumed.net

These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located. Specific questions physicians may have about the availability and use of the products described on these materials should be directed to their particular authorized Acumed distributor. Specific questions patients may have about the use of the products described in these materials or the appropriateness for their own conditions should be directed to their own physician.