


OSTEOMED

A TRADITION OF QUALITY AND INNOVATION

SPECTRUM™

LeFORT III

MID-FACE DISTRACTION SYSTEM



- THE INTERNAL DESIGN OF SPECTRUM™ LeFORT III DOES NOT REQUIRE AN EXTERNAL HALO DEVICE OR DISTRACTION WIRES RESULTING IN EASE OF SURGEON USE
- THE INTERNAL DESIGN PROVIDES ENHANCED PATIENT QUALITY OF LIFE AND SAFETY DURING THE DIFFERENT PHASES OF DISTRACTION
- ACTIVATION POINTS ARE EASILY ACCESSIBLE
- THE SPECTRUM™ LeFORT III CAN BE USED FOR EITHER A LeFORT III DISTRACTION OR A COMBINATION OF LeFORT I AND III DISTRACTIONS

PRODUCTS THAT IMPROVE PEOPLE'S LIVES

SPECTRUM™

Refer to *Product Insert* for specific indications and *Instructions For Use*.

LEFORT III MID-FACE DISTRACTION SYSTEM

VARIOUS
DISTRACTION
ROD LENGTHS

MULTIPLE CRANIAL
ATTACHMENT OPTIONS

FRONT PORTION CAN
BE USED FOR STABILITY
ON SIMULTANEOUS
LEFORT I DISTRACTION

PART NUMBER	DESCRIPTION	PART NUMBER	DESCRIPTION
217-0055	Spectrum™ Distraction Tool, Horizontal Adjustment	211-1604	1.6mm x 4mm Auto-Drive™ Screw
217-0057	Spectrum™ Distraction Wrench, Vertical Adjustment	211-1605	1.6mm x 5mm Auto-Drive™ Screw
217-0060	Spectrum™ Short Malar Anchor Assembly	211-1606	1.6mm x 6mm Auto-Drive™ Screw
217-0061	Spectrum™ Long Malar Anchor Assembly	211-1608	1.6mm x 8mm Auto-Drive™ Screw
217-0062	Spectrum™ Left Maxillary Plate Assembly	211-2004	2.0mm x 4mm Auto-Drive™ Screw
217-0063	Spectrum™ Right Maxillary Plate Assembly	211-2005	2.0mm x 5mm Auto-Drive™ Screw
217-0064	Spectrum™ Small Horizontal Bow	211-2006	2.0mm x 6mm Auto-Drive™ Screw
217-0065	Spectrum™ Large Horizontal Bow	211-2008	2.0mm x 8mm Auto-Drive™ Screw
217-0066	Spectrum™ Set Screw	220-0061	1.6mm Screwdriver Shaft, Manual
217-0067	Spectrum™ Distraction Rod Cap	220-0021	2.0/2.4mm Screwdriver Shaft, Manual
217-0068	Spectrum™ Small Straight Horizontal Bow		
217-0069	Spectrum™ Large Straight Horizontal Bow		
217 0081	Spectrum™ LeFort III, Anchor Plate Assembly, 3mm		
217 0082	Spectrum™ LeFort III, Anchor Plate Assembly, 5mm		
217 0083	Spectrum™ LeFort III, Anchor Plate Assembly, 8mm		
217 0084	Spectrum™ LeFort III, Malar Pin		
217 0085	Spectrum™ LeFort III, Short Activation Rod		
217 0086	Spectrum™ LeFort III, Medium Activation Rod		
217 0087	Spectrum™ LeFort III, Long Activation Rod		
217 0088	Spectrum™ LeFort III, Cap Nut		

* Optional pilot drills available for 1.6mm and 2.0mm screws.



3885 Arapaho Road
Addison, Texas 75001
Customer Service: 800.456.7779
Main: 972.677.4600
Fax: 800.390.2620

www.osteomedcorp.com

P/N 030-1517 Rev.A



A COLSON ASSOCIATE