

Acu-Loc 2 with the Avulsion Hook Plate Marginal Rim Fracture Fixation

The Acu-Loc 2 Wrist Plating System offers one of the most comprehensive plating solutions in distal radius and ulna fracture fixation on the market. When used with the sterile-packed Avulsion Hook Plate, surgeons can more easily address the challenges of marginal rim fractures in volar and dorsal compartments.

Avulsion Hook Plate

- ▶ Extended fixation and articular buttress
- ▶ Ultra-thin 0.8 mm profile
- ▶ Versatile design enables both volar and dorsal approaches for marginal rim fracture fixation



Volar (Alt.)



Dorsal (with Frag-Loc)



Learn more about Acu-Loc 2
and the Avulsion Hook Plate
go.acumed.net/AL2-AHP



Acu-Loc® Wrist Plating System
Trusted Over One Million Times

Comprehensive, Intraoperative Versatility

Acu-Loc 2 offers a wide array of plates and screw technology specifically designed to address simple to complex distal radius and ulna fractures.

Clinical Reliability

Trusted by over one million surgeons worldwide, Acu-Loc 2 is anatomically contoured for consistent plate placement and a robust, rigid construct.


Innovative Technology

Patented Frag-Loc[®] Compression Screw and reduction tools such as radiopaque targeting guides and kickstands assist in fragment fixation



Key Features

- ▶ Low, 1.6 mm profile at distal edge of the standard volar distal radius (VDR) plate helps minimize flexor pollicis longus (FPL) complications
- ▶ Two volar distal radius plating families provide surgical options in plate placement
- ▶ Two divergent radial styloid screws offer two points of fixation for stability as well as up to 30° condylar angulation for increased stability in targeting the fragment
- ▶ Optional Avulsion Hook Plate helps capture marginal rim fragments via volar or dorsal surgical approaches



Learn more about Acu-Loc 2
and the Avulsion Hook Plate
go.acumed.net/AL2-AHP

These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located. Specific questions physicians may have about the availability and use of the products described on these materials should be directed to their particular authorized Acumed distributor. Specific questions patients may have about the use of the products described in these materials or the appropriateness for their own conditions should be directed to their own physician.

