

# OsteoMatch™

PATIENT MATCHED IMPLANTS

- **Optimal anatomical fit versus conventional fixation**
- **Reduce operating time**
- **Excellent aesthetic results**
- **Autoclavable**
- **Fracture resistant**



*Rethinking Possibilities, Reshaping Lives*



# OsteoMatch™ Product Realization

OsteoMatch PEEK Patient-Specific Cranial Implants are individually sized and shaped implantable prosthetic cranioplasty plates intended to correct defects/replace voids in the cranial skeleton of a specific patient. OsteoMatch is digitally designed and manufactured from a patient CT.

The app delivers a new tool for doctors to remotely educate, plan, review and share with both patients and colleagues. All cases are stored on our cloud server for easy and secure access.



**\* Critical Path:** Must be completed by Dr. or rep to proceed