



*Rethinking Possibilities, Reshaping Lives*

510K  
CLEARED

THE  
**P** **INNACLE**  
OF PRECISION

# PINNACLE

The Pinnacle Battery Powered Smart Driver features an ergonomic streamlined pencil grip and an optimized center of gravity for supreme control and comfort. OsteoMed's torque-limiting software sets Pinnacle apart from the competition. Programmed based on an exhaustive screw characterization, the software senses when any OsteoMed screw is fully seated and stops to prevent stripping of the bone. With a top speed of 2,000 RPM, Pinnacle seats neuro screws in less than a second. This smart power driver is 510(k) cleared.




# PINNACLE



## Mode Selection

Hold down back mode button for 2 seconds to switch between drive and drill modes.



## Speed and Precision

With a top speed of 2,000 RPM, Pinnacle seats neuro screws in less than a second and is also powerful enough to drill pilot holes into dense bone.

## Pinnacle Accessories

Pinnacle accepts OsteoMed hex shank attachments including driver stems and twist drills.



## Forward / Reverse Buttons



## Pencil Grip

Streamlined pencil grip for supreme control and comfort.

## Disposable Batteries

Batteries come in single or double sterile packed boxes.



## Pinnacle Ordering Information

Pinnacle Driver: 110-1001

Single Sterile Battery Pack: 110-1000-SP

Double Sterile Battery Pack: 110-1000-SP-02

Pinnacle Single Sterilization Tray: 110-1011

Pinnacle Double Sterilization Tray: 110-1012


Hex Shank Driver Stems: 114-XXXX

Hex Shank Twist Drills: 111-XXXX

## OsteoMed Products




DuraKnife



OsteoMed Mesh



Fast-Flap




OsteoVation




OsteoDriver 2




Profile Zero



OsteoForm Mesh



SmartFlex



OsteoMatch



*Rethinking Possibilities, Reshaping Lives*

3885 Arapaho Rd.

Addison, TX 75001

Customer Service: 800.456.7779

Outside the U.S.: 001.972.677.4600

Fax: 800.390.2620

Fax Outside the U.S.: 001.972.677.4709

E-mail: [customer.service@osteomed.com](mailto:customer.service@osteomed.com)

[www.osteomed.com](http://www.osteomed.com)

P/N 030-2004 Rev.A