



**LOW PROFILE**  
MIDFACE SYSTEM



*Rethinking Possibilities, Reshaping Lives*

# LOW PROFILE MIDFACE SYSTEM

- Accepts 1.6mm Profile Zero screws, to create a low profile plate to screw construct
- Beveled edges to minimize palpability
- Meshes available in two profile heights- .5mm & .6mm - as well as rigid versions
- Plates easily contourable to patient mid-face anatomy

## STANDARD



- .8 Plate & Screw Thickness

## LOW PROFILE



- .5 Plate & Screw Thickness
- Beveled Edges
- Reduced Palpability

## PLATES



213-1603 Straight Plate with Gap, 4 Hole



213-1606 Straight Plate, 8 Hole



213-1611 Straight Plate, 24 Hole



213-1607 Square Plate, Large




213-1608 5x2 Hole, Large



213-1615 Square Plate, Small



213-1616 5x2 Hole, Small



213-1602 Y Plate



213-1600 Double Y




213-1609 Curved Plate, 4 Hole




213-1610 Curved Plate, 8 Hole



213-1604 L Plate, Left




213-1605 L Plate, Right



213-1601 T Plate

## MESH



- 213-1612 93x93x0.5
- 213-1613 93x93x0.6
- 213-1614 203x203x0.5

## DRILLS & DRIVERS



220-0091 1.3mm PILOT DRILL, J-LATCH 3.5mm STOP ASSEMBLY



220-0263 1.6mm DRIVER STEM, PROFILE 0



220-0209 1.3mm PILOT DRILL, J-LATCH



Rethinking Possibilities, Reshaping Lives

3885 Arapaho Road  
Addison, Texas 75001  
Customer Service: 800.456.7779  
Main: 972.677.4600  
Fax: 800.390.2620  
www.osteomedcorp.com

