



A New Option

in the Acu-Loc® Family



acumed®

Acu-Loc® Wrist Spanning Plate

+ acumed | Acu-Loc® Wrist Spanning Plate

Acumed now offers the Acu-Loc Wrist Spanning Plate for complex distal radius fractures. The plate is designed to hold the wrist in distraction and to provide ligamentotaxis to the wrist on a temporary basis while the distal radius heals.



LOCKING AND NONLOCKING HEXALOBE SCREWS

The Acu-Loc Wrist Spanning Plate uses 2.7 mm and 3.5 mm Locking and Nonlocking Hexalobe screws for plate fixation.

Plates are designed to be placed over the second or third metacarpal

Optional distal radius screw cluster for additional fixation

Precontoured with a 6° dorsal bend

Proximal slot allows up to 2 mm of additional distraction

Two different plate lengths to accommodate varying patient anatomy

Length: 171 mm

Metacarpal Holes
2.7 mm Locking and Nonlocking Holes

Distal Radius Holes
2.7 mm Locking Holes

Internal fixation achieved through three incision points

Radial Shaft Holes
3.5 mm Locking Holes

3.5 mm Nonlocking Distraction Slot

Length: 188 mm

For more information, please visit:
go.acumed.net/wsp-1



HNW40-06-A

Effective: 6/2014

© 2014 Acumed® LLC

Acumed® Headquarters
5885 NW Cornelius Pass Road
Hillsboro, OR 97124

Office: 888.627.9957

Fax: 503.520.9618

acumed.net

These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located. Specific questions physicians may have about the availability and use of the products described on these materials should be directed to their particular local sales representative. Specific questions patients may have about the use of the products described in these materials or the appropriateness for their own conditions should be directed to their own physician.

Acumed® and Acu-Loc® are registered trademarks of Acumed, LLC.